ScanX Swift View

EN-
US

ES

FR

Installation and Operating Instructions
Instrucciones para montaje y uso
Notice de montage et d'utilisation

RxOnly

2142100008L41

AIR
TECHNIQUES equipped for life®

G8505

2202V007

ScanX Swift View

EN-
US

Installation and Operating Instructions

RxOnly

2142100008L41

AIR
TECHNIQUES equipped for life®

G8505

2202V007

Contents

Important information

1	About this document	3
1.1	Warnings and symbols	3
1.2	Copyright information	4
2	Safety	5
2.1	Unauthorized modification	5
2.2	Intended use (FDA)	5
2.3	General safety information	5
2.4	Specialist personnel	5
2.5	Protection from electric shock	5
2.6	FCC note	6
2.7	Notification requirement of serious incidents	6
2.8	Only use genuine parts	6
2.9	Transport	6
2.10	Disposal	7
2.11	Protection from cybersecurity threats	7

Product description

3	Overview	8
3.1	Scope of delivery	9
3.2	Accessories	9
3.3	Optional items	9
3.4	Consumables	9
3.5	Wear parts and spare parts	9
4	Technical data	10
4.1	Phosphor storage plate scanner	10
4.2	Phosphor storage plate	12
4.3	Barrier envelope	12
4.4	Model identification plate	13
4.5	Conformity assessment	13
5	Function	13
5.1	Phosphor storage plate scanner	13
5.2	Phosphor storage plate	14
5.3	Barrier envelope	15
5.4	Stylus	15

5.5	Protective cover (optional)	15
-----	---------------------------------------	----

Installation

6	Requirements	16
6.1	Installation/setup room	16
6.2	System requirements	17
6.3	Monitor	17
7	Installation	18
7.1	Setting up the unit	18
7.2	Removing the protective film from the touch screen	18
7.3	Attaching the Stylus	18
7.4	Checking the memory card	18
7.5	Electrical connection	19
7.6	Connecting the unit to the network	19
8	Commissioning and first start-up	21
8.1	Configuring the network	21
8.2	Configuring the unit in VisionX	21
8.3	Configuring the unit in DBSWIN	22
8.4	X-ray unit settings	24
8.5	Commissioning tests	25

Usage

9	Operating the touch screen	26
9.1	Navigating	26
9.2	Using menus	26
9.3	Entering text in a field	27
9.4	Calling up messages on the touch screen	27
10	Correct use of phosphor storage plates	28
11	Operation	29
11.1	X-ray	29
11.2	Scanning the image data via a computer	31
11.3	Scanning image data via the touch screen on the unit	32

- 11.4 Erasing the imaging plate 34
- 11.5 Switching the unit off 34
- 12 Cleaning and disinfection 35**
 - 12.1 Phosphor storage plate scanner . 35
 - 12.2 Barrier envelope 36
 - 12.3 Phosphor storage plate 36
 - 12.4 Stylus 37
- 13 Maintenance 38**
 - 13.1 Recommended maintenance
schedule 38
 - 13.2 Check image quality 38

Troubleshooting

- 14 Tips for operators and service technicians 40**
 - 14.1 Poor X-ray image 40
 - 14.2 Software error 44
 - 14.3 Fault on the unit 45
 - 14.4 Error messages on the touch
screen 46

Appendix

- 15 Menu structure Settings 49**
- 16 Scanning times 52**
- 17 File sizes (uncompressed) 53**
- 18 Handover record 54**

Important information

1 About this document

These installation and operating instructions are an integral part of the unit.

Air Techniques shall not be held liable and offers no guarantees of the safe and smooth operation of this unit if you fail to comply with notes and instructions contained in these Installation and Operating Instructions.

The German version of the installation and operating instructions is the original manual. All other languages are translations of the original manual. These operating instructions are valid for ScanX Swift View, order number: G8500 (2142100510).

1.1 Warnings and symbols

Warnings

The warning notes in this document highlight possible injury to persons or damage to machinery.

They are marked with the following warning symbols:

General warning symbol

The warnings are structured as follows:

SIGNAL WORD

Description of type and source of danger

Here you will find the possible consequences of ignoring the warning

- Follow these measures to avoid the danger.

The signal word differentiates between different levels of danger:

- **DANGER**
Direct danger of severe injury or death
- **WARNING**
Possible danger of severe injury or death
- **CAUTION**
Risk of minor injuries
- **NOTICE**
Risk of extensive material/property damage

Miscellaneous symbols

These symbols are used in the document and on or in the unit:

Note, e.g. specific instructions regarding the efficient use of the unit.

Part number

Serial number

Medical device

Health Industry Bar Code (HIBC)

Batch name

CE mark

CE labeling with the number of the notified body

Manufacturer

Date of manufacture

Distributor

Dispose of correctly in accordance with EU Directive 2012/19/EU (WEEE).

Take note of the accompanying electronic documents.

Please note the operating instructions.

Wear hand protection.

Disconnect all power from the unit.

Do not reuse

DC current

-  Warning – risk of dangerous electric voltages
-  Warning - laser beams
-  This way up / store and transport in an upright position
-  Keep dry
-  Stacking limits
-  Lower and upper humidity limits
-  Lower and upper temperature limits
-  Upper temperature limit
-  Fragile, handle with care
-  Keep away from sunlight during storage

Adhesive label


Fig. 1: Laser class 3B


Fig. 2: Warning - laser beams


Fig. 3: Specification of laser source

1.2 Copyright information

All electronic drawings, processes, names, software, and appliances mentioned here are protected under copyright. Printing or copying these Installation and Operating Instructions, including excerpts thereof, may only be carried out with the written approval of Air Techniques.

2 Safety

The unit has been developed and designed appropriately such that hazards are largely excluded if the unit is used in accordance with its Intended use.

Therefore, please note the following.

Despite this, the following residual risks can remain:

- Personal injury due to incorrect use/misuse
- Personal injury due to mechanical effects
- Personal injury due to electric shock
- Personal injury due to radiation
- Personal injury due to fire
- Personal injury due to thermal effects on the skin
- Personal injury due to lack of hygiene, e.g. infection

2.1 Unauthorized modification

Pursuant to Part 15.21 of the FCC rules, any changes or modifications to this equipment not expressly approved by the manufacturer may cause, harmful interference and void the FCC authorization to operate this equipment.

2.2 Intended use (FDA)

ScanX Swift View

Indications for use

The ScanX Swift View is intended to be used for scanning and processing digital images exposed on Phosphor Storage Plates (PSPs) in dental applications.

Contraindications

Any other usage or usage beyond this scope is deemed to be improper. The manufacturer accepts no liability for damage resulting from improper usage. The user bears the sole risk. This unit is not suitable for monitoring patients over extended periods of time. This unit must not be used in operating rooms or similar rooms, in which hazards may arise from the combustion of flammable materials. The touch screen only shows a preview as an initial impression of the X-ray image. For purposes of diagnosis, the X-ray image must be viewed on a diagnostic monitor. The preview of the X-ray image on the touch screen is not suitable for the purposes of diagnosis.

Barrier envelope

Indication for use

Disposable Barrier Envelopes are intended to be used as a disposable barrier for Air Techniques Phosphor Storage Plates. This device is non-sterile and intended for single patient use only.

Contraindication

This Device has no contraindications.

2.3 General safety information

The sale or prescription of this device by a medical practitioner is subject to the restrictions of the applicable Federal Acts. The device may be used only under permanent supervision by a dentist or licensed medical practitioner.

Rx_{only} Caution: By virtue of Federal Law, the device may only be sold to dentists or bought on behalf of a dentist.

- › Comply with the guidelines, laws, rules and regulations applicable at the site of operation when you use this unit.
- › Prior to each use, check the function and proper condition of the device.
- › Do not convert or modify the unit.
- › Comply with the Installation and Operating Instructions.
- › Make the Installation and Operating Instructions always available to the operator in the vicinity of the device.

2.4 Specialist personnel

Operation

Persons operating the unit must ensure safe and correct handling based on their training and knowledge.

- › Instruct or have every user instructed in handling the unit.

Installation and repairs

- › All installation, resetting, alteration, expansion, and repair work must be carried out either by Air Techniques personnel or by a suitably qualified person approved by Air Techniques.

2.5 Protection from electric shock

- › Comply with all relevant electrical safety regulations when you work with this unit.

- › Never touch the patient and unshielded plug connections of the device at the same time.
- › Replace any damaged cables or plugs immediately.

Comply with the EMC rules concerning medical devices

The device is designed for use in a basic electromagnetic environment in accordance with IEC 61326-1.

NOTICE

Negative effects on the EMC due to non-authorized accessories

- › Use only Air Techniques accessories or accessories approved by Air Techniques.
- › Using any other accessories may result in increased electromagnetic interference emissions or the unit having reduced electromagnetic immunity, leading to an erroneous operation mode.

NOTICE

Erroneous operation mode due to use immediately adjacent to other devices or with other stacked devices

- › Do not stack the unit together with other devices.
- › If this is unavoidable, the unit and other devices should be monitored in order to ensure that they are working correctly.

2.6 FCC note

This device has been tested and complies with the limits for a class B digital device as set out in part 15 of the FCC regulations. These limits are designed to provide adequate protection against harmful interference in residential environments. This device generates, uses, and can emit high-frequency energy, and it can cause harmful interference to wireless communication systems if it is not installed and used in accordance with the instructions. However, this statement does not offer a guarantee that no interference will occur in a certain installation. If this device causes harmful interference to radio or TV reception (this can be tested by switching the device off and back on again), users are recommended to attempt to

eliminate the interference by taking one or more of the following measures:

- › Realign or reposition the reception aerial.
- › Increase the distance between the device and the receiver.
- › To meet FCC requirements on HF exposure compliance, a safety distance of at least 20 cm must be maintained between the device's internal aerial and all persons.
- › Connect the device to a power outlet on a different circuit to the one the receiver is connected to.
- › Consult a dealer or experienced radio/TV technician.

2.7 Notification requirement of serious incidents

The operator/patient has to report any serious incident related to the product to the manufacturer and the competent authority of the Member State, in which the operator and/or patient is established/resident.

2.8 Only use genuine parts

- › Only use accessories and optional items specified or approved by Air Techniques.
- › Only use original working parts and spare parts.

Air Techniques accepts no liability for damage resulting from the use of non-approved accessories, optional items or any parts other than original spare and wear parts.

The use of non-approved accessories, optional items or non-genuine wear parts / replacement parts (e. g. mains cable) can adversely affect the electrical safety and EMC.

2.9 Transport

Only the original packaging ensures optimum protection for the unit during transport. If necessary, the original packaging for this unit can be ordered from Air Techniques.

Air Techniques cannot be held responsible for any damage resulting from transport in unsuitable packaging, even during the warranty period.

- › Only transport the unit in its original packaging.

- › Keep all packaging away from children.
- › Do not expose the unit to any strong vibrations or shocks.

2.10 Disposal


Disposal of the units, electronic circuitry and PSPs must be accomplished only at the appropriate facilities for recovery and recycling. Make sure to dispose of such items in accordance with current federal, national, state and local government rules and regulations.

2.11 Protection from cybersecurity threats

The unit is to be connected to a computer that can be connected to the Internet. Therefore, the system needs to be protected from threats from the Internet.

- › Use antivirus software and update it regularly.
Look for evidence of possible virus infection and, if applicable, check with the antivirus software and remove the virus.
- › Perform regular data backups.
- › Provide access to units only to trustworthy users, e.g. by means of user name and password.
- › Make sure that only trustworthy contents are downloaded. Install manufacturer-authenticated software and firmware updates only.

3 Overview

- 1 ScanX Swift View PSP scanner
- 2 Stylus
- 3 Intraoral phosphor storage plate (PSP)
- 4 Barrier envelope intraoral
- 5 SDHC memory card (installed in the unit)
- 6 Network cable
- 7 WLAN stick
- 8 Power supply unit with country-specific adapter

3.1 Scope of delivery

The following items are included in the scope of delivery (possible variations due to country-specific requirements and/or import regulations):

ScanX Swift View

PSP scanner G8500A

- ScanX Swift View basic unit
- Power supply unit
- Network cable
- SDHC memory card
- Stylus
- VisionX Imaging software data carrier/download
- DBSWIN Imaging software data carrier/download
- ScanX View third-party software instructions data carrier/download
- WLAN stick (pre-installed depending on version)
- Phosphor storage plates:
 - Size 0
 - Size 2
- Barrier envelopes:
 - Size 0
 - Size 2
- PSP cleaning wipes
- Installation and Operating Instructions
- Quick Reference Guide

3.2 Accessories

The following items are required for operation of the device, depending on the application:

SDHC memory card G8590

Phosphor storage plates (PSPs)

Intraoral phosphor storage plate size 0
(2 pcs.) 73445-0

Intraoral phosphor storage plate size 1
(2 pcs.) 73445-1

Intraoral phosphor storage plate size 2
(4 pcs.) 73445-2

Intraoral phosphor storage plate size 3
(2 pcs.) 73445-3

Intraoral phosphor storage plate, size 4
(1 pcs.) 73445-4

Barrier envelopes

Barrier envelope Size 0
2 x 3 cm (100 pcs.) G8511-0

Barrier envelope Size 1
2 x 4 cm (100 pcs.) G8511-1

Barrier envelope Size 2
3 x 4 cm (300 pcs.) G8511-2

Barrier envelope Size 2
3 x 4 cm (1000 pcs.) G8511-2k

Barrier envelope Size 3
2.7 x 5.4 cm (100 pcs.) G8511-3

Barrier envelope Size 4
5.7 x 7.6 cm (50 pcs.) G8511-4

3.3 Optional items

The following optional items can be used with the device:

Network cable (5 m) G8595

Stylus G8585

Protective cover G8423

Commissioning and intraoral constancy tests

2D X-ray test phantom G8795

3.4 Consumables

The following materials are consumed during operation of the device and must be re-ordered:

Cleaning and disinfection

PSP cleaning wipes (50 pieces) B8910

Barrier envelopes

refer to "3.2 Accessories"

3.5 Wear parts and spare parts

Phosphor storage plates (PSPs)

refer to "3.2 Accessories"

4 Technical data

4.1 Phosphor storage plate scanner

Electrical data of the device

Voltage	V DC	24
Max. current consumption	A	1.25
Output	W	< 30
Type of protection		IP20

Electrical data of the power supply unit

Nominal Voltage	V AC	100 - 240
Frequency	Hz	50/60

General technical data

Dimensions (W x H x D)	mm	226 X 275 X 243
	in	8.90 X 10.83 X 9.57
Weight	kg	approx. 7
	lbs	approx. 15
Duty cycle S2 (in accordance with IEC 60034-1)	min	25
Duty cycle S6 (in accordance with IEC 60034-1)	%	25
Pixel size (selectable)	µm	12.5 - 50
Max. theoretical resolution	Line pairs/mm (Lp/mm)	approx. 40

Noise level

Ready to scan	dB(A)	approx. 37
During scanning	dB(A)	approx. 55

Network connection

LAN technology		Ethernet
Standard		IEEE 802.3u
Data rate	Mbit/s	100
Connector		RJ45
Type of connection		Auto MDI-X
Type of cable		≥ CAT5

WiFi connection

WiFi technology		IEEE 802.11b/g
Encryption		WPA, WPA2

Ambient conditions during operation

Temperature	°C	+10 to +35
	°F	+50 to +95
Relative humidity	%	20 - 80
Air pressure	hPa	750 - 1060

Ambient conditions during storage and transport

Temperature	°C	-20 to +60
	°F	-4 to +140
Relative humidity	%	10 - 95
Air pressure	hPa	750 - 1060

Classification

Medical Device Class (MDR)	I
FDA classification (CFR Title 21)	II
Laser class (unit) in accordance with IEC 60825-1:2014	1

Laser source

Laser class In accordance with IEC 60825-1:2014	3B	
Wavelength λ	nm	635
Output	mW	10

Memory card

Type	SDHC	
Maximum memory capacity	GB	32
File system	FAT32	
Performance class	Class	≥ 4

4.2 Phosphor storage plate

Classification

Medical Device Class (MDR)	Ila
FDA classification (CFR Title 21)	I

Ambient conditions during operation

Temperature	°C	18 - 45
	°F	64 - 113
Relative humidity	%	< 80

Ambient conditions during storage and transport

Temperature	°C	< 45
	°F	< 113
Relative humidity	%	< 80

Dimensions of intraoral phosphor storage plates

Size 0	mm	22 x 35
	in	0.87 x 1.38
Size 1	mm	24 x 40
	in	0.94 x 1.57
Size 2	mm	31 x 41
	in	1.22 x 1.61
Size 3	mm	27 x 54
	in	1.06 x 2.13
Size 4	mm	57 x 76
	in	2.24 x 2.99


4.3 Barrier envelope

Classification

Medical Device Class (MDR)	I
FDA classification (CFR Title 21) Class	II

4.4 Model identification plate

The type plate is located on the rear of the device.


REF Order number
SN Serial number

4.5 Conformity assessment

This device has been subjected to conformity acceptance testing in accordance with the current relevant guidelines of the European Union. This equipment conforms to all relevant requirements.

5 Function

5.1 Phosphor storage plate scanner

- 1 Operating elements
- 2 Input unit
- 3 Stylus
- 4 Release key
- 5 Collection tray
- 6 Memory card slot

The phosphor storage plate scanner is used to read image data stored on the phosphor storage plate.

The unit can be used in two different ways: via the imaging software (e. g. VisionX) on a PC or directly via the touch screen on the unit.


The transport mechanism guides the phosphor storage plate through the device. A laser in the scanner unit scans the PSP. The scanned data is converted into a digital image.

If a scanning job is started via the imaging software, the image is automatically transmitted to the computer.

If a scanning job is started via the touch screen, the image is saved to the memory card and needs to be transferred to the computer later on. After scanning, the phosphor storage plate runs through the erasure unit. Image data still held on the phosphor storage plate is erased with the aid of bright light.

The phosphor storage plate is then ejected for re-use.

Operating elements


- 1 Touch screen
- 2 On/off switch

The touch screen allows the unit to be operated when it is not connected to a computer. Instructions can be entered on the touch screen either with the tip of a finger or the stylus.

The **Help** button can be used to open a help page for the respective page. The **Messages** button can be used to call up current messages.

Connections

The connections are located on the rear of the unit, underneath the cover.

- 1 Connection for power supply unit
- 2 Reset button
- 3 USB connection (for WLAN stick)
- 4 Network connection status LEDs
- 5 Network connection
- 6 AUX connection for diagnostic units

ScanManager

When the ScanManager is enabled, more than one X-ray job can be transmitted simultaneously to the unit from different computers. The unit manages the X-ray jobs in a queue from which the respective X-ray job can be selected using the touch screen and then executed.

Without ScanManager, the unit processes one X-ray job at a time and is blocked until this job has been completed. During this time, no further X-ray jobs can be transmitted to the unit from other computers.

 ScanManager can be enabled via **Settings > System Settings > Operating Type**.


5.2 Phosphor storage plate

The PSP stores X-ray energy, which is re-emitted in the form of light after excitation by the laser.

This light is then converted into image information in the phosphor storage plate scanner.

The PSP has an active side and an inactive side. The PSP must always be exposed on the active side.

When used properly, a PSP can be exposed, read and erased several hundred times provided there is no mechanical damage. The PSP must be replaced if there are any signs of damage, e.g. if the protective layer is damaged or there are visible scratches that could interfere with the diagnostics.

- | | | |
|---|---------------|--|
| 1 | Inactive side | Black, size and manufacturer information printed on it |
| 2 | Active side | Light blue, with positioning aid a |

Positioning aid **a** is visible on the x-ray image and makes orientation easier during diagnosis.

5.3 Barrier envelope


Disposable Barrier Envelopes are intended to be used as a disposable barrier for Air Techniques Phosphor Storage Plates. This device is non-sterile and intended for single patient use only.

5.4 Stylus

The touch screen can be operated using the Stylus as an alternative to the tip of a finger.

5.5 Protective cover (optional)

The protective cover protects the device against dust and dirt, e. g. during extended periods, in which it is not in use.

Only qualified specialists or persons trained by Air Techniques may install, connect, and commission the unit.

6 Requirements

6.1 Installation/setup room

The room chosen for set up must meet the following requirements:

- Closed, dry, well-ventilated room
- It should not be a room made for another purpose (e.g. boiler room or wet cell)
- Max. light intensity 1000 Lux, no direct sunlight at the place of installation of the unit
- There should be no major fields of interference (e.g. strong magnetic fields) present that can interfere with the proper operation of the unit.
- Ambient conditions correspond to "4 Technical data".

6.2 System requirements

Interface:	Ethernet \geq 100 Mbit/s WLAN IEEE 802.11b/g with WPA and WPA2 encryption
Software:	VisionX version 2.3 or higher (order number: E7300) VisionX Connect version 3.0 or higher DBSWIN version 5.14 or higher (order number: E7200A), VistaEasy, Image Bridge

For the system requirements of the computer systems, visit the download area at www.airtechniques.com (document no. E7201).

6.3 Monitor

The monitor must meet the requirements for digital X-ray with a high light intensity and wide contrast range.

Strong ambient light, sunlight impinging directly onto the monitor and reflections can make it more difficult or even impossible to perform a diagnosis based on the X-ray images.

7 Installation

7.1 Setting up the unit

NOTICE

Damage to sensitive components of the unit due to shocks or vibrations

- › Do not expose the unit to any strong vibrations or shocks.
- › Do not move the unit during operation.

Portable and mobile HF communication devices can interfere with the effectiveness of electrical medical devices.

The unit can be set up as a tabletop unit or mounted on a wall using the wall bracket.


The load-bearing capacity of the table or wall must be suitable for the weight of the unit (see "4 Technical data").

Setting the unit on a table

To prevent errors when scanning the image data, install the unit so it is not exposed to vibrations.

- › Place the unit on a firm, horizontal surface.


7.2 Removing the protective film from the touch screen

- › Grasp one corner of the protective touch screen film and peel it off carefully.

7.3 Attaching the Stylus

- › The stylus is held on the unit by a magnet. Place the stylus in the indentation provided.

7.4 Checking the memory card


CAUTION

Loss of image data due to unexpected insertion or removal of the memory card

- › Insert or remove the memory card only when unit is switched off.

- › Check whether the memory card has been inserted correctly in the unit. If the memory is inserted incorrectly in the device, take it out again and re-insert it properly.

7.5 Electrical connection

Electrical safety when making connections

- › Connect the device to a correctly installed power outlet only.
- › Do not operate any other systems using the same multiple socket.
- › Make sure that none of the electrical cables leading to the unit are under any mechanical tension.
- › Before initial start-up verify that the mains supply voltage and the voltage stated on the type plate match (see also "4. Technical data").


Connecting the unit to the mains supply

 The unit has no main power switch. For this reason, it is important to set up the unit properly such that the plug can be easily accessed and unplugged if required.


Requirements:

- ✓ Properly installed power outlet close to the unit (observe the max. length of the power cord)
- ✓ Easily accessible power outlet
- ✓ Mains voltage matches the information shown on the type plate of the power supply unit


- › Attach the matching country-specific mains cable to the power supply unit.

- › Remove the cover from the rear of the device.

- › Plug the connecting plug of the power supply unit into the connection socket of the unit.
- › Secure the cable with a cable clip.

- › Plug the mains plug into the power outlet.
- › Replace the cover.

 The cover on the rear must be correctly fitted when the device is operated within the patient environment.

7.6 Connecting the unit to the network

Purpose of the network connection

The network connection is used to exchange information or control signals between the unit and a software installed on a computer, in order to, e. g.:

- Display parameters
- Select operating modes
- Indicate messages and error situations
- Change device settings
- Activate test functions
- Transmit data for archiving
- Provide documents concerning the devices

The unit can be connected to the network with a network cable or via WLAN.

 For information on connection via WLAN, see "Configuring WLAN on the unit".

Combining devices safely


- Safety and essential performance features are independent of the network. The device is designed appropriately for operation independent of a network. However, some of the functions are not available in this case.
- Faulty manual configuration can lead to significant network problems. The expert knowledge of a network administrator is required for configuration.
- The data connection utilizes part of the bandwidth of the network. Interactions with other medical devices cannot be completely excluded. Apply the IEC 80001-1 standard for risk assessment.
- The device is not suitable to be connected directly to the public internet.

Take care when connecting units together or to parts of other systems as there is always an element of risk (e.g. due to leakage currents).

- › Only connect units when there can be no question of danger to operator or to patient.
- › Only connect units when it is safe to do so and when there is no risk of damage or harm to the surroundings.
- › If it is not completely clear from the data sheet of the unit that such connections can be safely made or if you are in any doubt, always get a suitably qualified person (e.g. the relevant manufacturer) to verify that the setup is safe.
- › Only connect peripheral units (e.g., computer, monitor, printer) that conform at least to the requirements set out in IEC 60950-1.

Connecting the unit via the network cable

- › Remove the cover from the rear of the device.
- › Plug the enclosed network cable into the network socket of the device.

- › Replace the cover.

 The cover on the rear must be correctly fitted when the device is operated within the patient environment.

8 Commissioning and first start-up

NOTICE

Short circuit due to the build-up of condensation

- › Do not switch on the unit until it has warmed up to room temperature and is dry.

The unit supports the following imaging programs:

- VisionX manufactured for Air Techniques
- DBSWIN manufactured for Air Techniques
- Please inquire about third-party software

Always use the current version of the imaging program in the commissioning of the device. Check the version of the enclosed imaging program versus the versions available at www.airtechniques.com.

8.1 Configuring the network

- › Switch on the network devices (router, PC, switch).
- › Check if TCP port 2006 and UDP port 514 are enabled in the firewall; enable them, if necessary.
- If you are using the Windows firewall, you do not need to check the ports since you will be asked whether you want to enable them during the driver installation process.

Network configuration

Various options are available for network configuration:

- ✓ Automatic configuration via DHCP.
- ✓ Automatic configuration via Auto-IP for direct connection of device and computer.
- ✓ Manual configuration.
- › Configure the network settings of the device using the software or, if applicable, the touch screen.
- › Check the firewall and release the ports, if applicable.

Network protocols and ports

Port	Purpose	Service
45123 UDP, 45124 UDP	Device recognition and configuration	

Port	Purpose	Service
2006 TCP	Device data	
514 ¹⁾ UDP	Event protocol data	Syslog

- 1) The port can vary depending on the configuration.

Configuring WLAN on the unit

If the unit is to be operated via WLAN, the enclosed WLAN stick must be inserted into the USB connection on the back of the unit (see "Connections"). The connection to the unit then needs to be configured.

Proper functioning is only ensured if the WiFi stick supplied by Air Techniques is used.

In order to establish a secure WiFi connection, we recommend encrypting the WiFi network with WPA2.

The quality and transmission range of the WiFi connection can be reduced by environmental conditions (e.g. thick walls, other WiFi devices). Mind the signal strength when you select a suitable location for set-up.

Requirements:

- ✓ You need to be logged-in on the unit as Administrator or Service Technician (**Settings > Access Levels > Administrator/Service Technician**).
- › Check the WiFi settings with you Network Administrator.
- › Tap the following on the touch screen: **Settings > System Settings > Network**.
- › Under **Interface**, select the **WLAN** option and confirm with **OK**.
- › Configure the WLAN.
- › Confirm with **OK**.

When the unit is first connected to a computer, it applies the language and time settings of the computer.

8.2 Configuring the unit in VisionX

Configuration is performed directly in VisionX.

- ›  **Select the units.**

Entering a permanent IP address (recommended)


To reset the network settings, keep the unit reset button pressed for 15 - 20 seconds while switching on.

- › Deactivate *DHCP*.
- › Enter the IP address, subnet mask and gateway.
- › Click on *Apply*.
This saves the configuration.

Testing the device

You can scan an X-ray image to check if the unit is properly connected.

- › Select the *Test* tab.

- › Select the unit from the *Registered Devices* selection list.
- › Select the mode class.
- › Select the mode.
- › Click on *Scan Image*.
- › For the scanning of an image plate, see "11 Operation".

8.4 X-ray unit settings

A setting of 60 kV is preferred provided it can be set on the X-ray unit. The standard exposure values for F-speed film (e. g. Kodak Insight) can be used.

Intraoral X-ray units for an adult patient

The following table shows the standard values for the exposure time and the dose area product of a PSP for an adult patient.

	DC emitter, 7 mA Tube length 20 cm					
	without X-ray field limitation		X-ray field limitation 0.8 x 1.2 in ² (20 x 30 mm ²)		X-ray field limitation 1.2 x 1.6 in ² (30 x 40 mm ²)	
	60 kV	mGycm ²	60 kV	mGycm ²	60 kV	mGycm ²
Incisor	0.08 s	14.6	0.08 s	3.1	0.08 s	6.2
Premolar	0.12 s	21.9	0.12 s	4.6	0.12 s	9.3
Molar	0.17 s	31.1	0.17 s	6.6	0.17 s	13.2
Bite wings	0.18 s	32.9	0.18 s	7.0	0.18 s	14

	DC emitter, 6 mA Tube length 30 cm					
	without X-ray field limitation		X-ray field limitation 0.8 x 1.2 in ² (20 x 30 mm ²)		X-ray field limitation 1.2 x 1.6 in ² (30 x 40 mm ²)	
	70 kV	mGycm ²	70 kV	mGycm ²	70 kV	mGycm ²
Incisor	0.13 s	11.8	0.13 s	2.5	0.13 s	5.0
Premolar	0.18 s	16.4	0.18 s	3.4	0.18 s	6.9
Molar	0.25 s	22.8	0.25 s	4.8	0.25 s	9.6
Bite wings	0.27 s	24.6	0.27 s	5.2	0.27 s	10.4

› Check and adjust the specific X-ray unit in accordance with the standard values.

Intraoral X-ray units for a pediatric patient

The following table shows the standard values for the exposure time and the dose area product of a PSP for a pediatric patient.

Children are more sensitive to radiation than adults. Keep the exposure parameters as low as possible taking into consideration the image quality. Refer also to the FDA information on pediatric X-ray imaging.
(refer to the FDA website on this topic, see <https://tinyurl.com/FDAPediatric>).

	DC emitter, 7 mA Tube length 20 cm					
	without X-ray field limitation		X-ray field limitation 0.8 x 1.2 in ² (20 x 30 mm ²)		X-ray field limitation 1.2 x 1.6 in ² (30 x 40 mm ²)	
	60 kV	mGycm ²	60 kV	mGycm ²	60 kV	mGycm ²
Incisor	0.05 s	9.1	0.05 s	1.9	0.05 s	3.8
Premolar	0.07 s	12.8	0.07 s	2.7	0.07 s	5.4
Molar	0.11 s	20.1	0.11 s	4.2	0.11 s	8.5
Bite wings	0.11 s	20.1	0.11 s	4.2	0.11 s	8.5

	DC emitter, 6 mA Tube length 30 cm					
	without X-ray field limitation		X-ray field limitation 0.8 x 1.2 in ² (20 x 30 mm ²)		X-ray field limitation 1.2 x 1.6 in ² (30 x 40 mm ²)	
	70 kV	mGycm ²	70 kV	mGycm ²	70 kV	mGycm ²
Incisor	0.08 s	7.3	0.08 s	1.5	0.08 s	3.1
Premolar	0.11 s	10.0	0.11 s	2.1	0.11 s	4.2
Molar	0.14 s	12.8	0.14 s	2.7	0.14 s	5.4
Bite wings	0.14 s	12.8	0.14 s	2.7	0.14 s	5.4

- › Check and adjust the specific X-ray unit in accordance with the standard values.

8.5 Commissioning tests

The required tests (e.g. acceptance test) must be done in accordance with local rules and regulations.

- › Find out which tests are required.
- › Carry out testing in accordance with local rules and regulations.

Acceptance check

 The Intra / Extra Digital test phantom, and possibly the corresponding test phantom holder as well, is required for acceptance tests with the PSP and sensor as receivers.

- › Before commissioning the unit, the acceptance test of the X-ray system must be carried out in accordance with national regulations.
The constancy tests, which must be carried out at regular intervals by the surgery personnel, are based on the results of the acceptance test.

Electrical safety checks

- › Carry out the electrical safety check according to national law (e. g. in accordance with IEC 62353).
- › Document the results.
- › Carry out and document the instruction and handover for the unit.

 A sample handover report is included in the attachment.


9 Operating the touch screen

NOTICE

Damage to the touch screen due to incorrect handling

- › Only operate the touch screen using your fingertips or the Stylus.
- › Do not use sharp objects (e.g. ball-point pen) to operate the touch screen.
- › Protect the touch screen from water.

- › Operate the touch screen by tapping it with a fingertip or the Stylus to select a menu or field.

- › For further information about any window, tap on the *Help* field.

9.1 Navigating

If the contents of the window cannot be completely displayed on the touch screen, a scroll bar appears.

- › Tap  or  to move the displayed section of the window.

9.2 Using menus

The menus integrated in a window contain additional commands that can be selected.

- › To open the menu, tap .



Fig. 4: Example: expanded menu

- › Select a command.

9.3 Entering text in a field

- › If a field requires input, tap in the field. This opens the keyboard window.

The unit continues to operate.

-  Information Information for the operator. The unit continues to operate.

-  Fault-free operation

- › Tap on *Messages*. This displays the message. If there are several messages, the most current with the highest priority is displayed first.
- › For more information about the message, tap on *Help*.

9.4 Calling up messages on the touch screen

The *Messages* view shows a history of all previous messages. Messages are subdivided into the following categories:

-  Fault Unit will no longer function. When the error has been remedied, it may be necessary to acknowledge the error message.
-  Notice After acknowledgment the unit will continue to work, but only with limited functions.
-  Note Important information for the operator, e.g. about the current status of the device.

10 Correct use of phosphor storage plates

WARNING

Risk of cross contamination when not using the barrier envelope or when using the barrier envelope more than once

- › Do not use an phosphor storage plate without a barrier envelope.
- › Do not re-use the barrier envelope (disposable item).

CAUTION

Image data on the phosphor storage plate is not permanent

The image data is altered by light, natural X-ray radiation and scattered X-ray radiation. This will lead to a reduction in diagnostic information and clarity.

- › Read the image data within 30 minutes of exposure.
- › Never handle exposed phosphor storage plates without the barrier envelope.
- › Do not subject an exposed phosphor storage plate to X-ray radiation before and during the scanning process. Do not X-ray during the scanning process if the unit is in the same room as the X-ray tube.
- › Phosphor storage plates must only be read using a phosphor storage plate scanner that is approved by Air Techniques.


CAUTION

Phosphor storage plates (PSP) are toxic


PSPs that are not used with a barrier envelope can lead to poisoning when placed in the mouth or swallowed.

- › Only place PSPs in the patient's mouth in a barrier envelope.
- › Do not swallow the PSP or parts of it.
- › If the PSP or parts of it have been swallowed, consult a specialist doctor immediately and remove the PSP.
- › If the barrier envelope was damaged in the patient's mouth, rinse the mouth thoroughly with copious amounts of water. Do not swallow the water in the process.

- › Phosphor storage plates are flexible like X-ray film. However, the phosphor storage plates should not be kinked.

- › Do not scratch the phosphor storage plates. Do not subject the phosphor storage plates to pressure from hard or pointed objects.

- › Do not soil the phosphor storage plates.

- › Protect the phosphor storage plates from sunlight and ultraviolet light.
Store phosphor storage plates in a barrier envelope or intraoral/extraoral plate protector of proper size.
- › Phosphor storage plates will be pre-exposed on exposure to natural radiation and stray x-ray radiation. Protect erased or exposed phosphor storage plates from X-ray beams.
If the phosphor storage plate has been stored for longer than one week, erase the phosphor storage plate prior to use.
- › Do not store phosphor storage plates under hot or moist conditions. Note the ambient conditions (see "4 Technical data").
- › When used properly, phosphor storage plates can be exposed, read and erased several hundred times provided there is no mechanical damage.
Replace the phosphor storage plate if there are any signs of damage, e.g. protective layer is damaged or visible scratches, that could interfere with the diagnosis.
- › Phosphor storage plates that have a production or packaging defect will be replaced by Air Techniques in the same quantity.
- › Clean phosphor storage plates properly (see "12 Cleaning and disinfection").

11 Operation

CAUTION

Image data on the phosphor storage plate is not permanent

The image data is altered by light, natural X-ray radiation and scattered X-ray radiation. This will lead to a reduction in diagnostic information and clarity.

- › Read the image data within 30 minutes of exposure.
- › Never handle exposed phosphor storage plates without the barrier envelope.
- › Do not subject an exposed phosphor storage plate to X-ray radiation before and during the scanning process. Do not X-ray during the scanning process if the unit is in the same room as the X-ray tube.
- › Phosphor storage plates must only be read using a phosphor storage plate scanner that is approved by Air Techniques.

11.1 X-ray

The procedure is described using a size 2 PSP as an example.

Required accessories:

- PSP
- Barrier envelope of the same size as the PSP

WARNING

Risk of cross contamination when not using the barrier envelope or when using the barrier envelope more than once

- › Do not use an phosphor storage plate without a barrier envelope.
- › Do not re-use the barrier envelope (disposable item).

WARNING


Danger from the re-use of products intended for single use

Single-use article is damaged after use and cannot be reused.


- › Dispose of single-use articles after use.

Preparing the X-ray

- ✓ The PSP has been cleaned.
- ✓ The phosphor storage plate is not damaged.
- › If using it for the first time or if it has been stored for over a week: erase the phosphor storage plate (see "11.4 Erasing the imaging plate").
- › Completely slide the phosphor storage plate into the barrier envelope. The black (inactive) side of the PSP must be visible.

- › Pull off the adhesive strip and close the barrier envelope tightly by pressing together firmly.

- › In the case of phosphor storage plates Plus size 4, place the bite protector around the light protection cover with phosphor storage plate if necessary.

Taking the X-ray image

NOTICE


Damage to the phosphor storage plate (PSP) caused by a sharp-edged holding system

- › Only use holding systems that do not damage the barrier envelope or the PSP.
- › Do not use holding systems with sharp edges.

Wear hand protection.

- › Place the PSP in the barrier envelope into the patient's mouth. In doing this, make sure that the active side of the phosphor storage plate faces the X-ray tube.

- › Set the exposure time and setting values on the X-ray unit (see "8.4 X-ray unit settings").
- › Record an X-ray image. The image data must be scanned within 30 minutes.

Preparing for scanning

CAUTION

Light erases the image data on the phosphor storage plate

- › Never handle exposed phosphor storage plates without either a barrier envelope or an intraoral/extraoral plate protector.


Wear hand protection.

- › Remove the PSP with the barrier envelope from the patient's mouth.

**WARNING****Contamination of the unit**

- › Clean and disinfect the barrier envelope and protective gloves, e. g. wipe with a clean cellulose cloth.

- › In the event of heavy soiling, e. g. from blood, dry clean the barrier envelope and protective gloves, e. g. wipe with a clean cellulose cloth.
- › Disinfect the barrier envelope and protective gloves with a disinfection wipe.
Alternatively, use a spray disinfectant on a soft, lint-free cloth.

- › Allow the barrier envelope to fully dry.
- › Pull off the protective gloves, disinfect and clean the hands.

**NOTICE****Powder from the protective gloves on the PSP damages the unit during scanning**

- › Completely clean all traces of the protective glove powder from your hands before handling the PSP.

- › Tear off the barrier envelope.

**11.2 Scanning the image data via a computer****Starting the phosphor storage plate scanner and software**

The reading-out process is described using the VisionXimaging software.

For further information regarding the use of the imaging software, refer to the relevant manual.

- › Press the on / off switch  to switch on the device.
- › Switch on the computer and monitor.
- › Start VisionX.
- › Select the patient.
- › Select the corresponding image type in the menu bar.
- › Select the unit.
- › Set the acquisition mode.
Recording starts directly.
- › If ScanManager is enabled, select the X-ray job on the touch screen of the unit.

Result:

The touch screen will display an animated visual symbol requesting insertion of the phosphor storage plate.

Insert the phosphor storage plate only once the bar above the animation is green.


Fig. 5: Example of the animation prompting for insertion of the phosphor storage plate

Scanning the PSP


i In order not to mix up X-ray images, only scan the X-ray images from the selected patient.

- Place the barrier envelope with the PSP centrally and straight onto the input unit. The opened side of the barrier envelope faces down, the inactive side of the PSP faces the operator.

The feeder mechanism moves forward automatically and takes hold of the barrier envelope with PSP.

- Slide the PSP out of its barrier envelope downwards into the device until the PSP is automatically drawn in.


The barrier envelope is held in place by the feeder mechanism and is not transported into the unit.

Scanning progress is displayed on the touch screen. The image data is automatically transmitted to the imaging software. After it has been scanned, the PSP is erased and drops into the collection tray.

- Save the X-ray image.
- Remove the empty barrier envelope.

- Remove the PSP and prepare it for recording another X-ray image.

11.3 Scanning image data via the touch screen on the unit

Starting the phospor storage plate scanner

No PC connection is needed for scanning the image data via the touch screen. The image data is stored locally on the memory card. In order to transfer the image data to the imaging software, the unit must be connected to a computer. There are two options for scanning via the touch screen:

Scan:

Before scanning the image data, the patient data and exposure settings of the image are acquired and then saved with the image data.

If no patient data and exposure settings of the image are entered, the image is saved to a folder with date and time.

Quick scan:

The image data is saved to a folder with the date and time and no additional information.

Use **Help** on the touch screen for further information on operating the unit via the touch screen.

Prerequisites:

- ✓ Memory card (SDHC, max. 32 GB) is inserted in the slot of the unit.
- Press  to switch the unit on.

Start scanning:

- On the touch screen, tap on *Scan*.
- Enter the patient data.
- Select the image settings and scanning mode.

An animation on the touch screen prompts for insertion of the imaging plate.

-  Insert the imaging plate only once the bar above the animation is green.

Fig. 6: Example of the animation prompting for insertion of the imaging plate

Start quick scan:

- On the touch screen, tap on *Quick scan*.
- Select the scanning mode.

An animation on the touch screen prompts for insertion of the imaging plate.

-  Insert the imaging plate only once the bar above the animation is green.


Fig. 7: Example of the animation prompting for insertion of the imaging plate

Scanning the PSP


-  In order not to mix up X-ray images, only scan the X-ray images from the selected patient.

- Place the barrier envelope with the PSP centrally and straight onto the input unit. The opened side of the barrier envelope faces down, the inactive side of the PSP faces the operator.

The feeder mechanism moves forward automatically and takes hold of the barrier envelope with PSP.

- Slide the PSP out of its barrier envelope downwards into the device until the PSP is automatically drawn in.


The barrier envelope is held in place by the feeder mechanism and is not transported into the unit.

Scanning progress is displayed on the touch screen. The image data is automatically saved to the memory card.

-  The touch screen only shows a preview as an initial impression of the X-ray image. For purposes of diagnosis, the X-ray image must be viewed on a diagnostic monitor.

After it has been scanned, the PSP is erased and drops into the collection tray.

- EN-US
- › Remove the empty barrier envelope.
 - › Remove the PSP and prepare it for recording another X-ray image.

Transmitting image data to the computer

X-ray images created using the touch screen of the device are saved to a memory card. These X-ray images can be imported into an imaging software via a network connection.

- › Connect the unit to the network.
- › Start the imaging software.
- › Start the image import via the imaging software (further information can be found in the manual of the imaging software).
- › Save the image data.

The image data on the memory card is erased automatically as soon as the transfer has been successfully completed.

11.4 Erasing the imaging plate

The image data is erased automatically after scanning.

If you do not want the image data to be erased, this function can be disabled for the current scanning process by selecting *Disable erasing light* on the touch screen of the unit.

The special *ERASE* mode only activates the erasure unit of the imaging plate scanner. No image data is scanned.

The PSP needs to be erased using the special mode in the following cases:

- The first time the PSP is used, or if it is stored for more than one week.
- Due to an error, the image data on the PSP has not been erased (software error message).

Erasing the PSP via a computer

- › Select the special *ERASE* mode in the software.
- › Scan the PSP (see "11.2 Scanning the image data via a computer").

Erasing the PSP via the touch screen

- › On the touch screen, tap on *Quick scan*.
- › Select the scanning mode *ERASE*.
- › Scan the PSP (see "11.3 Scanning image data via the touch screen on the unit").

11.5 Switching the unit off

- › Press the on/off switch \odot for 3 seconds. As soon as the unit has shut down it switches off completely. The touch screen is off.

After you switch off the unit, wait 10 s before switching the unit on again.

Use of a protective cover

The protective cover protects the device against dirt and dust during extended periods in which it is not used.


WARNING

Danger of suffocation

- › Store the protective cover out of the reach of children.

- › Pull the protective cover over the device until it is completely covered. Make sure that the markings are at the front.

- › Store the protective cover in a safe place when it is not in use.

12 Cleaning and disinfection

Unless specified otherwise, Air Techniques recommends Monarch disinfectant wipes for cleaning and disinfection of the unit and accessories. For a full list of approved cleaning agents, please contact Air Techniques.

Do not use the following disinfectant wipes:

- CaviWipes 1
- Discide Ultra

NOTICE

The use of unsuitable agents and methods can damage the unit and accessories

- › Only use the disinfection and cleaning agents specified or approved by Air Techniques and the EPA.
- › Comply with the operating instructions of the disinfectants and cleaning agents.

Wear hand protection.

12.1 Phosphor storage plate scanner

Surface of the unit

The surface of the unit must be cleaned and disinfected if it is contaminated or soiled.

NOTICE

Liquid can cause damage to the unit

- › Do not spray the unit with cleaning agents or disinfectants.
- › Make sure that liquid penetrates into the unit.
- › Remove any soiling with a soft, lint-free cloth that has been dampened with cold tap water.
- › Disinfect the surfaces with a disinfectant wipe. Alternatively, use a spray disinfectant on a soft, lint-free cloth. Comply with the operating instructions of the disinfectant.

Input unit

The input unit must be cleaned and disinfected if there are indications of contamination or visible soiling.


NOTICE


Heat can damage plastic parts

› Do not treat any part of the unit in a thermal disinfector or steam sterilizer.


- › Tap  on the touch screen.
The fixing mechanism moves into the cleaning position.
- › Press the release key and remove the cover upwards at the same time.

- › Remove the fixing mechanism by moving it upwards.

- › Clean the cover, fixing mechanism and inside parts with a soft, lint-free cloth that has been dampened with cold tap water.

- › Disinfect the cover, fixing mechanism and inside parts with a disinfection wipe. Alternatively, use a spray disinfectant on a soft, lint-free cloth. Comply with the operating instructions of the disinfectant. The cover and fixing mechanism can also be disinfected in an immersion disinfection system.
- › Remount the fixing mechanism.
- › Remount the cover.
- › Touch **OK** on the touch screen. The fixing mechanism moves into the starting position.

12.2 Barrier envelope

- › Clean the Barrier Envelope after being removed from the patient's mouth with a disinfectant wipe, such as Air Techniques Monarch Surface Wipes.
- › Allow the Barrier Envelope to completely dry prior to ejecting the PSP.

12.3 Phosphor storage plate

Cleaning and disinfectant wipes are unsuitable for the cleaning of phosphor storage plates and may damage them. Only use a cleaning agent that is compatible with the materials: Air Techniques recommends the PSP Cleaning Wipes (see "3.4 Consumables"). Only this prod-

uct has been subjected to material compatibility testing by the manufacturer.

**NOTICE****Heat or humidity damage the phosphor storage plate**

- › Do not sterilize the phosphor storage plate with steam.
 - › Do not disinfect the phosphor storage plate by immersion.
 - › Only use cleaning agents that are compatible with the materials.
-
- › Soiling on both sides of the phosphor storage plate should be cleaned off with a soft, lint-free cloth before each use.
 - › Remove persistent or dried soiling with the phosphor storage plate cleaning cloth. Comply with the instructions for use of the cleaning cloth.
 - › Allow the phosphor storage plate to completely dry before using it.

12.4 Stylus

The Stylus can be cleaned in the same way as the unit (see "12.1 Phosphor storage plate scanner").

13 Maintenance

13.1 Recommended maintenance schedule

Do not attempt to carry out any maintenance work.

The internal components of the unit are accessible only after parts are removed with special tools and may be opened and subjected to maintenance by authorized service technicians of dealers only.

Contact your local authorized Air Techniques dealer for any maintenance. Non-compliance with these instructions can lead to equipment damage and loss of the warranty.

Prior to working on the unit or in case of danger, disconnect it from the mains.

The recommended maintenance intervals are based on using the device for 50 intraoral images per day on 220 working days per year.

Maintenance interval	Maintenance work
Annually	› Visually inspect the device.
	› Check the PSPs for scratches and replace them, if necessary.
	› Check the belt drives, transport belts and springs, and replace them if necessary.
	› Remove dust and dirt from accessible parts.
	› Carry out a system check.
Every 2 years	› Change the wiper blade set.
	› Change the roller fixtures.
	› Change the drive belt.

13.2 Check image quality

To assure the image quality, the unit must be subjected to maintenance (see "13.1 Recommended maintenance schedule") and regular cleaning and disinfection (see "12 Cleaning and disinfection") and the image quality of the PSP and X-ray system need to be checked.

Inspection interval	Work
Daily / before each use	› Clean the PSP, if necessary.
	› Check the PSP for scratches. If there are scratches on the surface, take a homogeneous test image (see "Check the PSP with a homogeneous test image"), replace the PSP, if applicable.
	› Keep an eye on the image quality with each X-ray image taken, also refer to "14.1 Poor X-ray image".
Monthly	› Take a homogeneous test image of the PSP (refer to "Check the PSP with a homogeneous test image"). If scratches or artifacts are visible in the image that may possibly have an adverse effect on the diagnostics, replace the PSP.
Every 3 months	› Check X-ray system.
	› Take an X-ray image with a 2D X-ray test phantom. Check the image for homogeneity, resolution, contrast and artifacts, refer to "Check X-ray system".

Check the PSP with a homogeneous test image

Scratches on the surface of the PSP may be visible in the X-ray image and may impair the ability to diagnose the X-ray images. If scratches are visible in the image that may impair the diagnostics, the PSP must be discarded.


Fig. 8: PSP with scratches

- › Place the PSP on a level surface at a distance of approximately 30 cm from the X-ray tube. Doing this, make sure that the active side of the PSP faces the X-ray tube.
- › Set the exposure time and setting values on the X-ray unit for a molar x-ray image (see "8.4 X-ray unit settings").
- › Scan the PSP with the PSP scanner at high resolution.
- › Check the image for homogeneity.
No scratches visible in the image: PSP can still be used.
Scratch is visible in the image: discard the PSP.

Check X-ray system

To check the X-ray system, take an X-ray image with the test phantom (refer to "3.3 Optional items"). This can be used to check the image produced with the X-ray system for homogeneity, resolution, contrast and artifacts.

- › Take an X-ray image with the test phantom. Comply with the instructions for use of the test phantom.
- › Read the PSP.
- › Check the image for homogeneity, resolution, contrast and artifacts.
- › If errors are visible in the image, contact a service technician.

? Troubleshooting




14 Tips for operators and service technicians


 Any repairs above and beyond routine maintenance may only be done by suitably qualified personnel or by one of our service technicians.




 Prior to working on the unit or in case of danger, disconnect it from the mains.


14.1 Poor X-ray image

Error	Possible cause	Remedy
X-ray image does not appear on the monitor after scanning	PSP inserted improperly and inactive side was scanned	› Scan the PSP again immediately, making sure you feed it in correctly in the process.
	Image data on the PSP has been erased, e.g. by ambient light	› Always scan the image data of the PSP as soon as possible.
	Fault in the unit	› Contact technician.
	No image data on PSP, PSP not exposed	› Expose the PSP.
	X-ray unit is faulty	› Contact technician.
X-ray image too dark	X-ray dose too high	› Check X-ray parameters.
	Incorrect brightness/contrast settings in the software	› Adjust the brightness of the X-ray image in the software.
X-ray image too bright	Exposed PSP has been exposed to ambient light	› Always scan the image data of the PSP as soon as possible.
	X-ray dose too low	› Check X-ray parameters.
	Incorrect brightness/contrast settings in the software	› Adjust the brightness of the X-ray image in the software.
X-ray image only shadowy	The X-ray dose on the PSP was insufficient	› Increase X-ray dose.
	Amplification (HV value) is set too low in the software	› Increase amplification (HV value).
	Unsuitable scanning mode selected	› Select a suitable scanning mode.
	The setting for the threshold value is too high	› Reduce the threshold value.

Error	Possible cause	Remedy
<p>Top or bottom bulge in the X-ray image</p> 	<p>PSP inserted off-center and at an angle</p>	<ul style="list-style-type: none"> › Check the error code on the touch screen. › Insert the PSP centrally and straight.
<p>X-ray image is mirror-inverted</p>	<p>PSP exposed on the wrong side.</p>	<ul style="list-style-type: none"> › Insert the PSP correctly in the barrier envelope. › Position the PSP correctly.
<p>Round shadow on the X-ray image</p> 	<p>Plus ID PSP (with marker) exposed on the wrong side</p>	<ul style="list-style-type: none"> › When taking an X-ray, make sure that the active side faces towards the X-ray tube.
<p>Ghosting or double exposure on X-ray image</p> 	<p>PSP exposed twice PSP not sufficiently erased</p>	<ul style="list-style-type: none"> › Only expose the PSP once. › Check the erasure unit for proper function. › Inform a service technician, if the problem persists.
<p>X-ray image mirrored in one corner</p> 	<p>PSP creased during X-ray exposure</p>	<ul style="list-style-type: none"> › Do not bend the PSP.

Error	Possible cause	Remedy
<p>Shadow on the X-ray image</p> 	<p>PSP removed from the barrier envelope before scanning</p>	<ul style="list-style-type: none"> › Do not handle PSPs without a barrier envelope. › Store the PSP in a barrier envelope.
<p>X-ray image cut off, part missing</p> 	<p>A metal part of the X-ray tube is in front of the X-ray beam</p> <hr/> <p>Faulty edge masking in imaging software</p>	<ul style="list-style-type: none"> › Recording an X-ray image, make sure there are no metal parts between the X-ray tube and the patient. › Check X-ray tube. › Deactivate edge masking.
<p>Software unable to combine the data to make a complete image</p>	<p>The X-ray dose on the PSP was insufficient</p> <hr/> <p>Amplification (HV value) is set too low in the software</p> <hr/> <p>Unsuitable scanning mode selected</p> <hr/> <p>The setting for the threshold value is too high</p>	<ul style="list-style-type: none"> › Increase X-ray dose. › Increase amplification (HV value). › Select a suitable scanning mode. › Reduce the threshold value.
<p>X-ray image has strips on image</p>	<p>PSP has been pre-exposed, e.g. by natural radiation or stray X-ray radiation</p> <hr/> <p>Parts of PSP exposed to light during handling</p> <hr/> <p>PSP dirty or scratched</p>	<ul style="list-style-type: none"> › If the PSP has been stored for more than one week, erase the PSP prior to use. › Do not expose exposed PSPs to bright light. › Scan image data within half an hour after the exposure. › Clean the PSP. › Replace scratched PSP.
<p>Bright stripes in the scanning window</p>	<p>Too much incident ambient light during the scanning process</p>	<ul style="list-style-type: none"> › Darken the room. › Turn the unit such that no light is directly incident on the input unit.

Error	Possible cause	Remedy
<p>Horizontal, grey lines in the X-ray image, extending beyond the left and right image edge</p> 	<p>Transport slipping</p>	<p>› Clean the transport mechanism, replace the transport belts if necessary.</p>
<p>X-ray image is stretched lengthwise with bright, horizontal stripes</p> 	<p>Wrong barrier envelope or PSP used</p>	<p>› Only use original accessories.</p>
<p>X-ray image split vertically into two halves</p> 	<p>Dirt in the laser slit (e.g. hair, dust)</p>	<p>› Clean the laser slit.</p>
<p>X-ray image with small bright spots or clouding</p>	<p>Micro scratches on the PSP</p>	<p>› Replace the PSP.</p>

Error	Possible cause	Remedy
Lamination of the PSP detaches at the edge 	Wrong retainer system used	<ul style="list-style-type: none"> › Only use original PSP and film retainer systems.
	PSP handled incorrectly.	<ul style="list-style-type: none"> › Use the PSP correctly. › Comply with the operating instructions of the PSP and film retainer systems.

14.2 Software error

Error	Possible cause	Remedy
"Too much ambient light"	Unit is exposed to too much light	<ul style="list-style-type: none"> › Darken the room. › Turn the unit such that no light can directly enter into the entry slot.
"Incorrect power supply unit"	Incorrect power supply unit connected	<ul style="list-style-type: none"> › Use the enclosed power supply unit.
"Overtemperature"	Laser or erasure unit too hot	<ul style="list-style-type: none"> › Switch the unit off and allow it to cool.
"Erasure unit fault"	LED defective	<ul style="list-style-type: none"> › Contact technician.
Imaging software fails to recognize the unit	Unit not switched on	<ul style="list-style-type: none"> › Switch the device on.
	Connecting cable between unit and computer not correctly connected	<ul style="list-style-type: none"> › Check the connecting cable.
	Computer does not detect any connection to the unit	<ul style="list-style-type: none"> › Check the connecting cable. › Check the network settings (IP address and subnet mask).
	Hardware error	<ul style="list-style-type: none"> › Contact technician.
	The IP address of the unit is being used by another unit	<ul style="list-style-type: none"> › Check the network settings (IP address and subnet mask) and assign a unique IP address to each unit. › Inform a service technician, if the problem persists.

Error	Possible cause	Remedy
The unit does not appear in the selection list in VistaConfig	Unit is connected downstream of a router	<ul style="list-style-type: none"> › Configure the IP address without an intermediate router on the unit. › Reconnect the router. › Manually enter the IP address in VistaConfig and register the unit.
	The IP address of the unit is being used by another unit	<ul style="list-style-type: none"> › Check the network settings (IP address and subnet mask) and assign a unique IP address to each unit. › Inform a service technician, if the problem persists.
The unit appears in the VistaConfig selection list, but connection is not possible	Subnet masks of the computer and the unit do not match	<ul style="list-style-type: none"> › Check subnet masks, adjust if necessary.
Error message "E2490"	The connection to the unit was interrupted while the software was still attempting to communicate with the unit	<ul style="list-style-type: none"> › Restore the connection to the unit. › Repeat the process.

14.3 Fault on the unit

Error	Possible cause	Remedy
Unit does not switch on	No mains voltage	<ul style="list-style-type: none"> › Check the mains cable and plug connection and replace if necessary. › Check the power supply unit. › If the green status LED does not light up, replace the power supply unit. › Check the mains fuse in the building.
	On/off switch is defective	<ul style="list-style-type: none"> › Contact technician.
Unit switches off again after a short time	Mains cable or power supply unit plug not inserted correctly	<ul style="list-style-type: none"> › Check the mains cable and plug connections.
	Hardware defect	<ul style="list-style-type: none"> › Contact technician.
	Mains supply voltage too low	<ul style="list-style-type: none"> › Check the mains supply voltage.
Unit is on, but there is no display on the touch screen	Touch screen initialization fault	<ul style="list-style-type: none"> › Switch the unit off and on again.
	Touch screen brightness set too dark	<ul style="list-style-type: none"> › Update firmware. › Increase the brightness of the touch screen.
	Touch screen defective	<ul style="list-style-type: none"> › Contact technician.

Error	Possible cause	Remedy
Loud operating noises after switching on lasting for more than 30 seconds	Radiation deflector defective	› Contact technician.
Unit not responding	The unit has not yet completed the startup procedure	› After switching the unit on, wait 20 - 30 seconds for the startup procedure to be completed.
	Unit is blocked by the firewall	› Enable the ports for the unit in the firewall settings.
Network connection has been disconnected	WLAN stick is not inserted	› Insert the WLAN stick into the device.
	Distance to WLAN router too far	› Set up the unit closer to the WLAN router.
	Walls between WLAN router and unit are too massive	› Set up the unit closer to the WLAN router.
	Another WLAN network is affecting the operation of the unit's WLAN network	› Change the frequency range of the WLAN network.
	Connecting cable between unit and computer not correctly connected	› Check the connecting cable.
	IP address of the device is used by another device	› Check the network settings (IP address and subnet mask) and assign a unique IP address to each unit. › Inform a service technician, if the problem persists.

14.4 Error messages on the touch screen

Error	Possible cause	Remedy
Error code 1008	Connection interrupted	› Update the firmware.
Error code 1010	Temperature of unit too high	› Allow the unit to cool down.
		› Contact technician.
Error code 1022	Subassembly not initialized	› Fault in software, update the software if required.
		› Contact technician.
Error code 1024	Internal data communication fault	› Switch the unit off and back on again.
		› Update the firmware.
		› Darken the room.
		› Turn the unit so that no light can fall directly into the insertion slot.

Error	Possible cause	Remedy
Error code 1026	Incorrect acquisition mode	<ul style="list-style-type: none"> › Select a different acquisition mode › Inform a service technician. › Update the firmware. › Reset the scanning modes to the factory settings via the unit interface or the Imaging Software.
Error code 1100	Permitted time for scan process exceeded	<ul style="list-style-type: none"> › Contact technician. › Check the belt drive. › Check for blockage, remove PSP from unit.
Error code 1104	Erasure unit fault	<ul style="list-style-type: none"> › Contact technician. › Replace the erasure unit.
Error code -1116	Drive feed blocked	<ul style="list-style-type: none"> › Remove the blockage. › Contact technician.
Error code -1117	Feed position error	<ul style="list-style-type: none"> › Contact technician. › Check the feed (ease of movement, drive belts).
Error code -1118	Input unit cover open	<ul style="list-style-type: none"> › Close the cover. › Deactivate cleaning mode.
Error code -1121	Input unit fixing mechanism missing	<ul style="list-style-type: none"> › Insert the fixing mechanism. › Deactivate cleaning mode.
Error code 1153	Unit fault	<ul style="list-style-type: none"> › Switch the unit off and on again. › Update firmware.
Error code 1154	Internal data communication fault	<ul style="list-style-type: none"> › Switch the unit off and on again. › Update firmware.
Error code 1160	Final radiation deflector rotation speed not reached	<ul style="list-style-type: none"> › Contact technician. › Update firmware. › Replace the radiation deflector subassembly if the problem occurs regularly.
Error code -1171	Fault on laser	<ul style="list-style-type: none"> › Return the unit for repair.
Error code 1172	SOL sensor timeout Fault on the laser, SOL sensor or pentaprism assembly	<ul style="list-style-type: none"> › Inform a service technician. › Update the firmware.
Error code 10000	Unit is exposed to too much light	<ul style="list-style-type: none"> › Darken the room. › Turn the unit such that no light can be directly incident in the entry slot.
Error code 10009	Internal communication error warning; unit remains ready for operation	<ul style="list-style-type: none"> › Update the firmware.

Error	Possible cause	Remedy
Error code -10015	PSP fed in off-center	› Insert the PSP centrally.
Error code 10017	Unit shuts down	› Wait until the unit has shut down completely
Error code 2	System error during startup of the unit	› Switch the unit off and back on again. › Update the firmware.
Error code 78	Memory card full	› Transmit image data to the computer. › Insert an empty memory card.
	Fault during memory cleanup	› Press and hold the reset button while switching on the unit. › Update firmware. › Press and hold the reset button while switching on the unit.
Firmware not running	A firmware update has been carried out	› Switch the unit off and on again.
	Internal communication fault	› Switch the unit off and on again.
Settings (e.g. language) reset after unit restart	Faulty configuration file	› Update the firmware. › Reset the configuration to the factory settings and reconfigure.
Warning message during shutdown of the unit	Not a malfunction	› Update the firmware.

 Appendix

15 Menu structure Settings

Device information ¹	Device data
	Dealer information
	Report
Access levels ¹	User
	Administrator
	Service technician
	Factory technician

System settings ²	Language	German (DE) English (EN) ...	
	Date & time	Date Time	
	Network	MAC address	
		Name	
		Interface	LAN WLAN
		DHCP	
		IP address	
		Subnet mask	
		Gateway	
	Acquisition settings	Patient ID	
Last name			
First name			
Date of birth			
Male/female			
Pregnancy			
Comment			
X-ray station			
X-ray parameters			
Acquisition type	INTRA Child ...		
	X-ray stations	Room 1 Room 2 ...	
		Touch screen	Brightness Touch screen calibration
Device settings			Standby Standby time Menu fade-out time Autom. rotation Query radiation dose Patient ID scheme Service information Service interval
	Operating mode	ScanManager	

Service menu ³	Testing			
	Scanning modes	Edit scanning modes		
		Display scanning modes		
	Maintenance			
	Messages			
	Diagnosis	Statistics	Display statistics counter	
			Display error counter	
		Manipulate	Transport	
			Insert	
			Erasure unit	
			Pentaprism	
			PMT	
		Request sensor values	Sensors	
			Temperatures	
			Internal device voltages	
	Oscilloscope			
	Check touch screen			
Display test images				
Factory settings	Reset scanning modes			

- 1 visible from access level **Operator** or higher
- 2 visible from access level **Administrator** or higher
- 3 visible from access level **Service Technician** or higher

16 Scanning times

The scanning time is the time required for complete scanning of image data and depends on PSP format and pixel size.

The time to image depends mainly on the computer system used and its work load. Times stated are approximate.

Theoretical resolution (LP/mm)	24	17	10
Pixel size (μm)	21	29.4	50.5
Resolution (dpi)	1210	864	503
Intra Size 0	20 s	14 s	8 s
Intra Size 1	21 s	15 s	9 s
Intra Size 2	22 s	16 s	10 s
Intra Size 3	25 s	18 s	11 s
Intra Size 4	34 s	25 s	16 s

17 File sizes (uncompressed)

The file sizes depend on the PSP format and the pixel size. File sizes stated are approximate and have been rounded upwards.

Suitable compression methods can considerably reduce the file size without loss of data.

Theor. resolution (LP/mm)	24	17	10
Pixel size (µm)	21	29.4	50.5
Resolution (DPI)	1210	864	503
Intra Size 0	3.49 MB	1.78 MB	0.60 MB
Intra Size 1	4.35 MB	2.22 MB	0.75 MB
Intra Size 2	5.76 MB	2.94 MB	1.00 MB
Intra Size 3	6.61 MB	3.37 MB	1.14 MB
Intra Size 4	19.65 MB	10.02 MB	3.40 MB

18 Handover record

This document confirms that a qualified handover of the medical device has taken place and that appropriate instructions have been provided for it. This must be carried out by a qualified adviser for the medical device, who will instruct you in the proper handling and operation of the medical device.

Product name	Order number (REF)	Serial number (SN)

- Visual inspection of the packaging for any damage
- Unpacking the medical device and checking for damage
- Confirmation of the completeness of the delivery
- Instruction in the proper handling and operation of the medical device based on the operating instructions

Notes:

Name of person receiving instruction: **Signature:**

Name and address of the qualified adviser for the medical device:

Date of handover: **Signature of the qualified adviser for the medical device:**

--	--

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Air Techniques, Inc.