
# ProVecta 3D Prime

EN-  
US

FR

ES

Operating instructions  
Notice d'utilisation  
Instrucciones de uso

**Rx**Only

**AIR**  
**TECHNIQUES** equipped for life®

2210200496L41

A7766


1911V006

# ProVecta 3D Prime

EN-  
US

Operating instructions

**Rx**Only

**AIR**  
**TECHNIQUES** equipped for life®

2210200496L41

A7766

1911V006

# Contents

## Important information

<b>1</b>	<b>About this document</b>	3
1.1	Warnings and symbols	3
1.2	Copyright information	4
<b>2</b>	<b>Safety</b>	5
2.1	Intended use and indication	5
2.2	Normal use	5
2.3	Improper use	5
2.4	Contraindication	5
2.5	General safety information	5
2.6	Radiation protection	5
2.7	Specialist personnel	6
2.8	Protection from electric shock	6
2.9	Essential performance characteristics	6
2.10	Notification requirement of serious incidents	6
2.11	Only use genuine parts	7
2.12	Transport	7
2.13	Disposal	7
2.14	Protection from cybersecurity threats	7

## Product description

<b>3</b>	<b>Overview</b>	8
3.1	Scope of delivery	9
3.2	Optional items	9
3.3	Consumables	9
<b>4</b>	<b>Technical data</b>	10
4.1	X-ray tube performance data	11
4.2	Dimensions	16
4.3	Model identification plate	17
4.4	Conformity assessment	17
<b>5</b>	<b>Function</b>	18
5.1	Functional description	18
5.2	Operating elements	18
5.3	Status LED	19
5.4	Positioning aids	19

5.5	Exposure button	20
5.6	Memory card slot	20
5.7	Sensor window	20

## Usage

<b>6</b>	<b>Operating the touch screen</b>	21
6.1	Navigating	21
6.2	Using menus	21
6.3	Calling up messages on the touch screen	21
<b>7</b>	<b>Operation</b>	22
7.1	Switch the device on	22
7.2	Making settings in the imaging software	22
7.3	Inserting the positioning aids	29
7.4	Positioning the patient	33
7.5	Start a test run	36
7.6	Taking the X-ray image	37
7.7	Emergency stop switch	38
<b>8</b>	<b>Cleaning and disinfection</b>	40
8.1	Surface of the unit	40
8.2	Positioning aids	40
<b>9</b>	<b>Reprocessing</b>	42
9.1	Risk analysis and classification	42
9.2	Reprocessing procedure in accordance with EN ISO 17664	42
9.3	General information	43
9.4	Preparation at the operating location	43
9.5	Manual cleaning, intermediate rinsing, disinfection, drying	43
9.6	Automatic cleaning, intermediate rinsing, disinfection, final rinse, drying	44
9.7	Check for function	44
9.8	Steam sterilization	44
9.9	Issue clearance for the parts for sterilization	45
9.10	Storing parts for sterilization	45
<b>10</b>	<b>Maintenance</b>	46


# ProVecta 3D Prime

EN-  
US

FR

ES

Operating instructions  
Notice d'utilisation  
Instrucciones de uso

**Rx**Only

**AIR**  
**TECHNIQUES** equipped for life®

2210200496L41

 **Important information****1 About this document**

These installation and operating instructions are an integral part of the unit.

Air Techniques shall not be held liable and offers no guarantees of the safe and smooth operation of this unit if you fail to comply with notes and instructions contained in these Installation and Operating Instructions.

The German version of the operating instructions is the original manual. All other languages are translations of the original manual.

These operating instructions are valid for the Pro-Vecta 3D Prime, order number: A7761.

Refer to the separate installation instructions for information about assembly, installation and configuration of the unit (document no. A7767).

**1.1 Warnings and symbols****Warnings**

The warning notes in this document highlight possible injury to persons or damage to machinery.

They are marked with the following warning symbols:

General warning symbol

The warnings are structured as follows:

**SIGNAL WORD****Description of type and source of danger**

Here you will find the possible consequences of ignoring the warning

- › Follow these measures to avoid the danger.

The signal word differentiates between different levels of danger:

- **DANGER**  
Direct danger of severe injury or death
- **WARNING**  
Possible danger of severe injury or death
- **CAUTION**  
Risk of minor injuries
- **NOTICE**  
Risk of extensive material/property damage

**Miscellaneous symbols**

These symbols are used in the document and on or in the unit:

Note, e.g. specific instructions regarding the efficient use of the unit.

Part number

Serial number

Medical device

CE labeling with the number of the notified body

Manufacturer

Caution: Federal law restricts this device to sale by or on the order of a Doctor.

Dispose of correctly in accordance with EU Directive 2012/19/EU (WEEE).

Type BF applied part

Do not reuse

Not sterile


Sterilize with steam at 134 °C

Protective ground connection

Equipotential bonding

-  Fragile, handle with care
-  Lower and upper atmospheric pressure limits
-  Lower and upper temperature limits
-  Lower and upper humidity limits
-  Stacking limits
-  Recycling
-  Keep dry
-  This way up / store and transport in an upright position
-  Keep away from sunlight during storage
-  Refer to Operating Instructions.
-  Wear hand protection.
-  Wear eye protection.
-  Use a mask.
-  Wear protective clothing.
-  Disconnect all power from the unit.
-  Notice
-  Emergency stop switch
-  Warning – X-rays

 Warning – risk of dangerous electric voltages

 Warning – X-rays

## 1.2 Copyright information

All electronic drawings, processes, names, software, and appliances mentioned here are protected under copyright.

Printing or copying these Installation and Operating Instructions, including excerpts thereof, may only be carried out with the written approval of Air Techniques.

## 2 Safety

The unit has been developed and designed appropriately such that hazards are largely excluded if the unit is used in accordance with its Normal Use.

Despite this, the following residual risks can remain:

- Personal injury due to incorrect use/misuse
- Personal injury due to mechanical effects
- Personal injury due to electric shock
- Personal injury due to radiation
- Personal injury due to fire
- Personal injury due to thermal effects to skin
- Personal injury due to lack of hygiene, e.g. infection

### 2.1 Intended use and indication

ProVecta 3D Prime is computed tomography x-ray unit intended to generate 3D and panoramic X-ray images in dental radiography for adult and pediatric patients. It provides diagnostic details of the maxillofacial areas for a dental treatment.

The device is operated and used by physicians, dentists, and x-ray technicians.

Not intended for mammography use.

### 2.2 Normal use

The unit must only be used by dentists or dental assistants who have been trained in the use of X-ray radiation in accordance with the applicable legal requirements.

### 2.3 Improper use

Any use of this appliance/these appliances above and beyond that described in the Installation and Operating Instructions is deemed to be incorrect usage. The manufacturer cannot be held liable for any damage resulting from incorrect usage. The operator will be held liable and bears all risks.

### 2.4 Contraindication

The radiobiological effects of X-ray beams in tissue result in the following contraindications:

- Pregnancy
- Pre-existing illnesses preventing the recording of a CBCT image
- Absence of a justified indication

Exceptions can be formulated at the discretion of the physician.

## 2.5 General safety information

- › Comply with the guidelines, laws, rules and regulations applicable at the site of operation when you use this unit.
- › Prior to each use, check the function and proper condition of the device.
- › Do not convert or modify the unit.
- › Comply with the Installation and Operating Instructions.
- › Make the Installation and Operating Instructions always available to the operator in the vicinity of the device.

The sale or prescription of this device by a medical practitioner is subject to the restrictions of the applicable Federal Acts. The device may be used only under permanent supervision by a dentist or licensed medical practitioner.

**Rx<sub>only</sub>** Caution: Federal law restricts this device to sale by or on the order of a Doctor.

## 2.6 Radiation protection

- › Comply with all applicable X-ray protection regulations and X-ray protection measures.
- › Use the prescribed X-ray protection equipment.
- › In order to reduce the level of X-ray exposure, we recommend the use of bismuth, lead shielding or protective aprons, especially for children and teenagers.
- › The persons operating the equipment must keep away from the X-ray unit while the exposure is being taken. The minimum distance required by the law must be maintained.
- › Children and pregnant women must consult a doctor before recording an X-ray image.
- › No person other than the patient is permitted to be present in the radiation room without X-ray protection measures. In exceptional circumstances another person may be present to provide assistance, but this must not be a member of the surgery staff. When the exposure is being taken, make sure that you maintain visual contact with the patient and the unit and keep talking to the patient.
- › The radiation room must be lockable to prevent entry by unauthorized persons.
- › If a fault occurs, abort the exposure immediately by releasing the trigger button.

## 2.7 Specialist personnel

### Operation

Persons operating the unit must ensure safe and correct handling based on their training and knowledge.

- › Instruct or have every user instructed in handling the unit.

### Installation and repairs

- › All installation, resetting, alteration, expansion, and repair work must be carried out either by Air Techniques personnel or by a suitably qualified person approved by Air Techniques.

## 2.8 Protection from electric shock

- › Comply with all relevant electrical safety regulations when you work with this unit.
- › Never touch the patient and unshielded plug connections of the device at the same time.
- › Replace any damaged cables or plugs immediately.

### Comply with the EMC rules concerning medical devices

The unit meets the requirements according to IEC 60601-1-2:2014.

- › The unit is intended for use in professional healthcare facilities (in accordance with IEC 60601-1-2). If the unit is operated in any other environment, potential effects on the electromagnetic compatibility must be taken into account.
- › Do not use the device near HF surgical devices and MRI equipment.
- › Keep a minimum distance of 30 cm between the device and other electrical devices.
- › Note that cable lengths and cable extensions have effects on electromagnetic compatibility.
- › No maintenance measures are required to maintain the basic EMC safety.
- › The emissions characteristics of this device render it suitable for use in industrial environments and hospitals (CISPR 11, Class A). When used in a residential environment (which normally requires Class B in accordance with CISPR 11), this device may not provide adequate protection from radio communication services. The operator may need to take corrective measures such as relocating or reorienting the device.

### NOTICE

#### Negative effects on the EMC due to non-authorized accessories

- › Use only Air Techniques accessories or accessories approved by Air Techniques.
- › Using any other accessories may result in increased electromagnetic interference emissions or the unit having reduced electromagnetic immunity, leading to an erroneous operation mode.

### NOTICE

#### Erroneous operation mode due to use immediately adjacent to other devices or with other stacked devices

- › Do not stack the unit together with other devices.
- › If this is unavoidable, the unit and other devices should be monitored in order to ensure that they are working correctly.

### NOTICE

#### Reduced performance features due to insufficient distance between unit and mobile HF communication devices

- › Keep at least 30 cm distance between the unit (including parts and cables of the unit) and mobile HF communication devices (wireless units) (including their accessories such as antenna cables and external antennas).

## 2.9 Essential performance characteristics

The unit does not have any essential performance characteristics as set out in IEC 60601-1 section 4.3.

## 2.10 Notification requirement of serious incidents

The operator/patient has to report any serious incident related the product to the manufacturer and the competent authority of the Member State, in which the operator and/or patient is established/resident.

## 2.11 Only use genuine parts

- › Only Air Techniques accessories and special accessories or those approved by Air Techniques may be used.
- › Only use original spare and replacement parts.

Air Techniques accepts no liability for damage resulting from the use of non-approved accessories, special accessories or any parts other than original spare and wear parts.

The use of non-approved accessories, special accessories or non-genuine working parts / spare parts (e.g. mains cable) can have a negative effect on the electrical safety and EMC.

The following accessories may affect EMC:

- Mains cable (3.6 m; order no.: A7782)
- Exposure switch (order no.: A7801)

## 2.12 Transport

Only the original packaging ensures optimum protection for the unit during transport.

If necessary, the original packaging for this unit can be ordered from Air Techniques.

Air Techniques cannot be held responsible for any damage resulting from transport in unsuitable packaging, even during the warranty period.

- › Only transport the unit in its original packaging.
- › Keep all packaging away from children.
- › Reattach the transport locking devices.
- › Do not expose the unit to any strong vibrations or shocks.
- › Do not bump or pull the unit.

## 2.13 Disposal

### Device

The units and electronic circuits may be disposed of only in proper facilities for recycling and recovery. Make sure to dispose of these objects in compliance with the applicable laws and regulations of the federal, national, state and municipal authorities.

### X-ray emitter


The X-ray emitter contains a tube that is potentially capable of imploding, lead cladding and mineral oil.

## 2.14 Protection from cybersecurity threats

The unit is to be connected to a computer that can be connected to the Internet. Therefore, the system needs to be protected from threats from the Internet.

- › Use antivirus software and update it regularly. Look for evidence of possible virus infection and, if applicable, check with the antivirus software and remove the virus.
- › Perform regular data backups.
- › Provide access to units only to trustworthy users, e.g. by means of user name and password.
- › Make sure that only trustworthy contents are downloaded. Install manufacturer-authenticated software and firmware updates only.

### 3 Overview

- |   |   |    |   |
|---|---|----|---|
| 1 | 3D and panoramic X-ray unit                     | 8  | Head supports with cushion                |
| 2 | Wall bracket                                    | 9  | Hygienic protective covers for bite block |
| 3 | Frame grabber card                              | 10 | Bite block                                |
| 4 | USB dongle                                      | 11 | Adapter for bite block                    |
| 5 | USB stick with device-specific calibration data | 12 | Chin holder for maxillary joint image     |
| 6 | Small parts set (e. g. screws, nuts)            | 13 | Chin holder for edentulous jaws           |
| 7 | Exposure switch (with holder)                   | 14 | Chin holder for sinus image               |
|   |   | 15 | Fiber optic cable                         |
|   |   | 16 | Mains cable for permanent connection      |

### 3.1 Scope of delivery

The following items are included in the scope of delivery (possible variations due to country-specific requirements and/or import regulations):

- ProVecta 3D Prime . . . . . A7750
- Fiber optic cable 10 m / 33 ft
  - Exposure switch and holder
  - Holder for bite block
  - Bite blocks (3x)
  - Chin support for edentulous jaws
  - Chin holder for maxillary joint image
  - Chin holder for sinus image
  - Head supports with cushion
  - Wall bracket top (16 in)
  - Hygienic protective covers for bite block (100 pieces)
  - Small parts set (e. g. screws, nuts)
  - Various housing parts
  - Operating instructions
  - Installation instructions
  - PCI Express frame grabber card
  - USB dongle
  - USB stick containing device-specific calibration data

If the mains cable of this unit is damaged it must only be replaced by an original mains cable from the manufacturer.

### 3.2 Optional items

The following optional items can be used with the device:

- Bottom wall bracket set, 16 inch, short (wall/wall installation) . . . . . A7754
- Fiber optic cable 5 m . . . . . A7757
- Fiber optic cable 20 m . . . . . A7758
- Silicone pads for head support Plus . . . . . A7803
- Stand . . . . . A7815

#### Acceptance and consistency check

- 2D X-ray Pan phantom set . . . . . A7556
- Consistency check test phantom 3D set . . . . . A7759
- Acceptance test phantom 3D set . . . . . A7760
- Ball phantom . . . . . A7330
- 2D X-ray test phantom holder . . . . . A7366

### 3.3 Consumables

The following materials are consumed during operation of the device and must be re-ordered:

Hygienic protective cover bite block (100 pieces) . . . . . A7395

Silicone pads for head support Plus . . . . . A7803

#### Cleaning and disinfection

Monarch surface disinfection wipes (hundred pcs. - 7" x 9") . . . . . H6171

Monarch enzymatic cleaner (84.5 oz. bottle) . . . . . H6201

Information on replacement part is available at the portal for authorized specialist dealers:  
[www.airtechniques.com](http://www.airtechniques.com).

## 4 Technical data

### Electrical data for the device

Nominal voltage	V AC	200 - 240
Frequency	Hz	50/60
Protection class		I
Operating mode X-ray tube		S6 = 6.3% 320 s duty cycle 20 s / 5 min (switch-on/switch-off time)
Operating mode height adjustment		S3 = 9% duty cycle 1 min / 9 min (switch-on/switch-off time)
Rated power	W	170
Maximum power	kVA	2.2
Fuses* (2 pcs.)	A	T 10.0 AH / 250 V~ (IEC 60127-2, Sheet 5)

\* The fuses may only be replaced by Air Techniques or by a body authorized by Air Techniques.

### Classification

Medical device of class	IIb
FDA classification (CFR Title 21)	II

### General technical data

Dimensions (W x D)	mm	572.5 x 1181 ± 12
	in	22.54 x 46.50 ± 0.47
Height	mm	1406 - 2206
	in	55.35 x 86.85
Weight	kg	180
	lbs	397

### Ambient conditions during storage and transport

Temperature	°C	-10 to +60
	°F	14 to 140
Relative humidity	%	10 - 75
Air pressure	hPa	860 - 1060

### Ambient conditions during operation

Temperature	°C	10 - 35
	°F	50 - 95
Relative humidity	%	30 - 75
Air pressure	hPa	860 - 1060

### X-ray emitter

Model		DG-07C11T2
Rated power	kW	1.6 (for 1 s)
Type of high-voltage generator		Inverter

<b>X-ray emitter</b>		
Nominal voltage, high-voltage generator	kV	50 - 99 ( $\pm 10\%$ )
Nominal current, high voltage generator	mA	4 - 16 ( $\pm 10\%$ , max. 75 kV 16 mA, max. 99 kV 10 mA)
Cooling, high-voltage generator		Automatic monitoring Shut-off at $\geq 60^\circ\text{C}$
Additional filtering at 50 kV	mm Al in Al	2.0 + 3.0 (added automatically for CBCT) 0.08 + 0.12 (added automatically for CBCT)
Integrated filtering at 50 kV	mm Al in Al	0.8 0.03
Total filtering at 50 kV	mm Al in Al	2.8 + 3.0 (added automatically for CBCT) 0.11 + 0.12 (added automatically for CBCT)
X-ray tube model		Toshiba D-052SB
Focal spot size as per IEC 60336 X-ray tube	mm in	0.5 0.02
Anode angle	$^\circ$	5
Anode heat capacity	kJ	35
Pulse/break ratio		1:60 or more
Duration of radiation exposure	s	0.5 - 20
Maximum current-time product per hour	mA s	960 (at 75 kV/16 mA)


#### 4.1 X-ray tube performance data

- Maximum deviation of the peak voltage from the displayed value  $\pm 10\%$
- Maximum deviation of the tube current from the displayed value  $\pm 20\%$
- Maximum deviation of the exposure time from the displayed value  $\pm 10\%$
- The device complies with the standards IEC 60601-1, IEC 60601-1-3 and IEC 60601-2-63.
- The lowest possible load factor is obtained with a combination of the settings 60 kV and 4 mA.

## Maximum Rating Charts


DC (Center Grounded)

Constant potential high-voltage generator  
Nominal Focus Spot Value: 0.5


## Emission and Filament Characteristics

Constant potential high-voltage generator  
Nominal Focus Spot Value: 0.5

## Anode Thermal Characteristics

**Detector Pano/CBCT**

Brand	Xmaru1404CF	
Type	CMOS photodiode array	
Pixel size	$\mu\text{m}$	49.5 99 (2x2 binning) 198 (4x4 binning)
		0.001949 0.003898 (2x2 binning) 0.007795 (4x4 binning)
Sensor size	mm	230 x 160 x 26
	in	9.06 x 6.30 x 1.02
Active surface area	mm	135.8 x 36.4
	in	5.35 x 1.43
Frame rate	fps	53.5
		107 (2x2 binning)
		308 (4x4 binning)
Greyscale	bit	14

Acquisition mode	FDD mm in	FOD mm in	ODD mm in	Image capture scale (magnification factor)
Panoramic	600 23.62	477.7 18.81	122.3 4.81	1.26

FDD: distance from focal spot to detector

FOD: distance from focal spot to object

ODD: distance from object to detector ( $\text{ODD} = \text{FDD} - \text{FOD}$ )

Image capture scale =  $\text{FDD}/\text{FOD}$

**Electromagnetic compatibility (EMC)  
Interference emission measurements**

High-frequency emissions in accordance with CISPR 11	Group 1
Interference voltage at the power supply connection CISPR 11:2009+A1:2010	Class A
Electromagnetic interference radiation CISPR 11:2009+A1:2010	Class A
Emission of harmonics IEC 61000-3-2:2005+A1:2008+A2:2009	Class A
Voltage changes, voltage fluctuations, and emission of flicker IEC 61000-3-3:2013	Conforms

**Electromagnetic compatibility (EMC)  
Interference immunity measurements on cover**

Immunity to interference by discharge of static electricity IEC 61000-4-2:2008 ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Conforms
Immunity to interference by high-frequency electromagnetic fields IEC 61000-4-3:2006+A1:2007+A2:2010 3 V/m 80 MHz - 2.7 GHz 80 % AM at 1 kHz	Conforms
Immunity to interference by near fields of wireless HF communication devices IEC 61000-4-3:2006+A1:2007+A2:2010 See Table on immunity to interference table by near fields of wireless HF communication devices.	Conforms

**Immunity to interference levels by near fields of wireless HF communication devices**

Radio service	Frequency band MHz	Test level V/m
TETRA 400	380 - 390	27
GMRS 460 FRS 460	430 - 470	28
LTE band 13, 17	704 - 787	9
GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE band 5	800 - 960	28

**Immunity to interference levels by near fields of wireless HF communication devices**

Radio service	Frequency band MHz	Test level V/m
GSM 1800 CDMA 1900 GSM 1900 DECT LTE bands 1, 3, 4, 25 UMTS	1700 - 1900	28
Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE band 7	2400 - 2570	28
WLAN 802.11 a/n	5100 - 5800	9

**Electromagnetic compatibility (EMC)****Interference immunity measurements on supply input**

Immunity to interference by rapid transient bursts – AC

voltage grid

IEC 61000-4-4:2012

Conforms

± 2 kV

100 kHz repetition frequency

Immunity to interference, line-line

IEC 61000-4-5:2005

Conforms

± 0.5 kV, ± 1 kV

Immunity to interference by surges, line-earth

IEC 61000-4-5:2005

Conforms

± 0.5 kV, ± 1 kV ± 2 kV

Immunity to interference, line-conducted disturbances

induced by high-frequency fields – AC voltage grid

IEC 61000-4-6:2013

3 V

0.15 - 80 MHz

Conforms

6 V

ISM frequency bands

0.15 - 80 MHz

80 % AM at 1 kHz

Immunity to interference due to voltage dips, short inter-  
ruptions and voltage fluctuations



IEC 61000-4-11:2004

0%  $U_T$  for 0.5 cycle


Conforms

0%  $U_T$  for 1 cycle70%  $U_T$  for 25/30 cycles (50/60 Hz)0%  $U_T$  for 250/300 cycles (50/60 Hz)

## 4.2 Dimensions

### 4.3 Model identification plate


- 1 Type plate of the unit
- 2 Type plate of the X-ray tube

### 4.4 Conformity assessment

This device has been subjected to conformity acceptance testing in accordance with the current relevant guidelines of the European Union. This equipment conforms to all relevant requirements.

## 5 Function


- 1 Emergency stop switch
- 2 On/off switch
- 3 X-ray tube
- 4 C-shaped elbow
- 5 Status LED
- 6 Operating elements
- 7 Memory card slot

### 5.1 Functional description

Similar to computer tomography or magnetic resonance tomography, sectional images can be generated with CBCT. With CBCT, an X-ray tube and an imaging sensor opposite it rotate around a seated or standing patient. The X-ray tube rotates through 180°-540° and emits a conical X-ray beam. The X-rays pass through the region under investigation and are measured for image generation by a detector as an attenuated grey scale X-ray image. Here, a large series of two-dimensional individual images is acquired during the revolution of the X-ray tube. Using a mathematical calculation on the rotating image series via a reconstruction computer, a grey value coordinate image is generated in the three spatial dimensions. This three-dimensional coordinate model corresponds to a volume graphic that is made up of individual voxels. This volume can be used to generate sectional images (tomograms) in all spatial dimensions as well as 3D views.

## 5.2 Operating elements

- 1 Touch screen
- 2 Button for opening/closing the head supports
- 3 Button for positioning beams on/off
- 4 Buttons for height adjustment

The touch screen can be used to operate the unit. Information can be entered on the touch screen with the tip of a finger.



Fig. 1: Monitor, unit ready to acquire image

- 1 Logged-in patient
- 2 Selected X-ray image
- 3 Display of the X-ray parameters (duration, DAP value, voltage and current)
- 4 Selected parameters

The **Help** button can be used to open a help page for the respective page. The **Messages** button can be used to call up current messages.


### 5.3 Status LED

The status LED uses different colors to display the different operating modes:

- Blue: unit ready for operation
- Green: unit ready to acquire image
- Yellow: X-ray beam active

### 5.4 Positioning aids

- 1 Lever for positioning the Frankfurt plane positioning beam
- 2 Frankfurt plane of the X-ray positioning beam
- 3 Head supports with cushion
- 4 Positioning aid, e. g. chin support with bite block
- 5 Upper canine positioning beam

- 6 Lever for positioning the upper canine positioning beam
- 7 Mid-sagittal positioning beam
- 8 Grips


The applied parts in accordance with IEC 60601-1 are:

- Grips
- Head supports with cushion
- Positioning aids (e.g. bite block and mounting for bite block, chin support for edentulous patients)


#### Description of the positioning aids

The positioning aids are used to correctly position the patient in the unit. The appropriate positioning aid is selected according to the selected image. The head supports gently keep the head of the patient in place.

#### Panoramic


Bite block and adapter bite block


Chin holder for edentulous jaws

Chin holder for maxillary joint image

Chin holder for sinus image


Head supports with cushion

## 5.5 Exposure button

### Exposure switch


The exposure switch is used to trigger the prepared image acquisition and start the X-ray exposure. The LED indicates the status of the unit, as does the LED on the unit.

- Green: unit ready for recording
- Yellow: X-ray beam active

- 1 Indicator lamp (LED)
- 2 Exposure button


## 5.6 Memory card slot

The unit has a slot for a memory card. The slot is needed for service purposes only.

## 5.7 Sensor window

The active sensor surface is indicated by the markers in the corners of the sensor window. The cross indicates the geometric mid-point of the active sensor surface.

- 1 Active sensor surface
- 2 Geometric mid-point of the active sensor surface

## 6 Operating the touch screen


### NOTICE

#### Damage to the touch screen due to incorrect handling

- › Only touch the touch screen with your fingertips.
- › Do not use sharp objects (e.g. ball-point pen) to operate the touch screen.
- › Protect the touch screen from water.

- › Operate the touch screen by tapping it with a fingertip to select a button or input field.

- › For further information about any window, tap on the *Help* field.

### 6.1 Navigating

If the contents of the window cannot be completely displayed on the touch screen, a scroll bar appears.

- › Tap  or  to move the displayed section of the window.

### 6.2 Using menus

The menus integrated in a window contain additional commands that can be selected.

- › To open the menu, tap .


Fig. 2: Example: expanded menu

- › Select a command.

### 6.3 Calling up messages on the touch screen

The *Messages* view shows a history of all previous messages. Messages are subdivided into the following categories:

#### Fault

Unit will no longer function. When the error has been remedied, it may be necessary to acknowledge the error message.

#### Notice

After acknowledgement the unit will continue to work, but only with limited functions.

#### Note

Important information for the operator, e.g. about the current status of the device. The unit continues to operate.

#### Information

Information for the operator. The unit continues to operate.

#### Fault-free operation

- › Tap on *Messages*. This displays the message. If there are several messages, the most current with the highest priority is displayed first.
- › For more information about the message, tap on *Help*.

## 7 Operation

### CAUTION

#### Health risks for the patient due to contraindications

- › Before using the unit on the patient, make sure that none of the listed contraindications are evident.

In this section, the term "Child" is used for children and adolescents from the age of 7.

### 7.1 Switch the device on


### CAUTION

#### Danger of injury by moving rotating unit arm


After the unit is switched on and the parameters are confirmed on the touch screen, the C-shaped angle connector piece is positioned. Persons can be injured during this.

- › Nobody may be present in the area of the rotating unit arm while the unit is being switched on.

- › Use the main power switch to switch the unit on.

The main power switch lights up green after it is switched on. The status LED lights up blue.

### 7.2 Making settings in the imaging software

The settings are described using the VisionX imaging software for an example. VisionX can be downloaded from [www.airtechniques.com](http://www.airtechniques.com).

For further information regarding the use of the imaging software, refer to the relevant manual.



Children are more sensitive to radiation than adults. Keep the exposure parameters as low as possible taking into consideration the image quality. Refer also to the FDA information on pediatric X-ray imaging.

(refer to the FDA website on this topic, see <https://tinyurl.com/FDAPediatric>).

### Preparing an X-ray image in VisionX

Requirements:

- ✓ VisionX is started.
- ✓ Patient is logged-in.
- ✓ No other image acquisitions are in progress (X-ray or video).
- › Click on the required image in the menu bar (e.g.  for a CBCT image).  
You can use  to call up further acquisition types that belong to the grouping e. g.  for 5x5 Maxilla Molar right (see "Image acquisition programs").

Depending on how the image acquisition types are configured, the acquisition of the X-ray image will start either immediately or you will first need to select an X-ray station.

- › If image acquisition does not begin automatically, select the X-ray station.  
The parameters, imaging volume and patient shape, are preselected according to the patient.
- › Check the parameters.  
Single-click on  to open the fly-out for setting the parameters. The changed parameters are immediately synchronized with the device.
- › If the preselected parameters are correct, continue to work directly on the unit.

### Parameter overview

Depending on the chosen image acquisition program, different parameters are available ( e. g. the image volume is not available for panoramic images, but the jaw arch is instead).

#### Image volume

The selected image volume influences the height of the image volume. The image volume "Child" has a reduced height. The diameter is identical.

Image volume "Normal"  
Size (W x H): approx. 100 x 85 mm

Image volume "Child"  
Size (W x H): approx. 100 x 70 mm

### Image quality

HQ image  
An improved signal/noise ratio is achieved via an extended exposure time.

SQ image  
This setting is used for standard images.

### Patient type

The selection of patient type depends on the patient's size and/or head circumference. This means that the preset patient type may need to be changed, if necessary.

The X-ray parameters are preset based on the patient type (see "Appendix").

If a child is selected, then the x-ray parameters are different:

- Reduced dose
- Shorter revolution time
- Smaller radiation field

Large patient

Average patient

Small patient

Child (< 13 years of age)

### Jaw arch

The selected jaw shape affects the rotational behavior of the rotating unit arm during image acquisition. This enables an image with an ideal layer position to be captured even on a particularly narrow or wide jaw.

Normal jaw arch

Narrow jaw arch

Wide jaw arch

Child's jaw arch

## Image acquisition programs

## CBCT images

**CBCT**

The CBCT image shows the jaw area.

The size of the jaw area shown depends on the selected image volume.

Resolution: 200  $\mu\text{m}$

R

**CBCT 5x5 Maxilla Molar right**

The X-ray image depicts the right molar region of the maxilla with a volume of 5x5 cm.

Resolution\*: 120  $\mu\text{m}$

R

**CBCT 5x5 Maxilla Premolar right**

The X-ray image depicts the right premolar region of the maxilla with a volume of 5x5 cm.

Resolution\*: 120  $\mu\text{m}$

**CBCT 5x5 Maxilla Front**

The X-ray image depicts the front region of the maxilla with a volume of 5x5 cm.

Resolution\*: 120  $\mu\text{m}$

L

**CBCT 5x5 Maxilla Premolar left**

The X-ray image depicts the left premolar region of the maxilla with a volume of 5x5 cm.

Resolution\*: 120  $\mu\text{m}$

L

**CBCT 5x5 Maxilla Molar left**

The X-ray image depicts the left molar region of the maxilla with a volume of 5x5 cm.

Resolution\*: 120  $\mu\text{m}$

R

**CBCT 5x5 Mandible Molar right**

The X-ray image depicts the right molar region of the mandible with a volume of 5x5 cm.

Resolution\*: 120  $\mu\text{m}$

R

**CBCT 5x5 Mandible Premolar right**

The X-ray image depicts the right premolar region of the mandible with a volume of 5x5 cm.

Resolution\*: 120  $\mu\text{m}$

**CBCT 5x5 Mandible Front**

The X-ray image depicts the front region of the mandible with a volume of 5x5 cm.

Resolution\*: 120  $\mu\text{m}$

L

**CBCT 5x5 Mandible Premolar left**

The X-ray image depicts the left premolar region of the mandible with a volume of 5x5 cm.

Resolution\*: 120  $\mu\text{m}$

L

**CBCT 5x5 Mandible Molar left**

The X-ray image depicts the left molar region of the mandible with a volume of 5x5 cm.

Resolution\*: 120  $\mu\text{m}$

\* The resolution can be changed to 80  $\mu\text{m}$  in the service menu of the unit.

For panoramic images of children, the size of the radiation field is reduced with the aid of an additional collimator. The radiation dose is significantly reduced for this image.

Panoramic images


**Standard**

The standard panoramic image records the complete dental area with ascending dental branches and maxillary joints.


**Front**

The image shows a reduced dental area without ascending dental branches.

**Right**

The image only shows the right dental area.


**Left**

The image only shows the left dental area.


**Orthogonal**

The image shows the complete dental area and is generated perpendicular to the maxillary arch. This prevents overlapping crowns.


**Bite wing**

The image shows the lateral dental area with a size limited to the bite wings.


**Bite wing front**

The image shows the anterior dental area with a size limited to the bite wings.

**Panoramic images****Bite wing right**

The image shows the right posterior region with a size limited to the bite wings.

**Bite wing left**


The image shows the left posterior region with a size limited to the bite wings.

**Temporomandibular imaging****Maxillary joint, lateral**

The image shows the lateral maxillary joints for an open and closed mouth in 4-fold depiction in one image.

**Maxillary joint, PA**

The image shows the posterior-anterior maxillary joints for an open and closed mouth in 4-fold depiction in one image.

**Sinus images****Sinus, lateral**

The X-ray image shows the lateral sinuses.

Sinus images

**Sinus, PA**

The X-ray image shows the posterior-anterior sinuses.

### 7.3 Inserting the positioning aids

For the acquisition of the X-ray image, the patient is positioned in the unit using the corresponding positioning aids and then accurately aligned using the positioning beams.

#### WARNING

#### Danger due to non-processed products

- Risk of infection for operator and patient
- › Reprocess the product correctly and sterilize it as required prior to the first use and after every subsequent use.
- › Do not reprocess disposable items.

#### WARNING


#### Danger from the re-use of products intended for single use

- Single-use article is damaged after use and cannot be reused.
- › Dispose of single-use articles after use.

#### Inserting the head supports

If no head supports are inserted or if they are dirty, insert new head supports before positioning the patient.


- › Remove any dirty head supports by pulling them out.
  - › Insert new head supports.
- When doing this, make sure that the cushions of the head supports face inwards.


#### Inserting the cushions of the head supports

If no cushions are inserted in the head supports or if they are dirty, insert new cushions before you position the patient.

- › Remove any dirty cushions by pulling them out.
- › Inserting the cushion holder.

- › Insert new cushions in the cutout provided on the head supports.


**EN-US** Inserting the positioning aid for CBCT images

We recommend using the holder for bite block on CBCT images. The bite block can be used optionally in addition to this.

In the case of edentulous patients, the chin support for edentulous jaws can be used.

The other positioning aids can also be used, depending on the application scenario.

- › Inserting the holder for bite block.

The bite block can be used with or without a hygienic protective cover.

We recommend using the bite block with a hygienic protective cover.

If the bite block is used without a hygienic protective cover, refer to the instructions under "Inserting the positioning aid for panoramic images with hygienic protective cover (optional)" and the reprocessing instructions under "9 Reprocessing".

**WARNING**

**There is a danger of cross contamination if hygienic protective covers are not used or if they are re-used**

- › Reprocess the bite block after using it without the hygienic protective cover.
- › Do not re-use the hygienic protective cover (disposable item).


**WARNING**

**Danger from the re-use of products intended for single use**


Single-use article is damaged after use and cannot be reused.

- › Dispose of single-use articles after use.

- › Optionally insert the bite block.

- › In the case of edentulous patients, the chin support for edentulous jaws can be used.

### Inserting the positioning aid for panoramic images

We recommend using the mounting for the bite block and the bite block for panoramic images. In the case of edentulous patients, the chin support for edentulous patients can be used. The other positioning aids can be used depending on the application scenario.

The bite block can be used with or without a hygienic protective cover.

We recommend using the bite block with a hygienic protective cover.


If the bite block is used without a hygienic protective cover, refer to the instructions under "Inserting the positioning aid for panoramic images with hygienic protective cover (optional)" and the reprocessing instructions under "9 Reprocessing".


### WARNING

**There is a danger of cross contamination if hygienic protective covers are not used or if they are re-used**

- › Reprocess the bite block after using it without the hygienic protective cover.
  - › Do not re-use the hygienic protective cover (disposable item).
- › Inserting the mounting for the bite block and the bite block.

- › In the case of edentulous patients, the chin support for edentulous patients can be used.

Inserting the positioning aid for panoramic images with hygienic protective cover (optional)


**WARNING**

**Risk of cross contamination due to non-reprocessed bite block**

› Reprocess the bite block in accordance with the reprocessing instructions.


› Optionally place a hygienic protective cover over the bite block.

**Inserting the positioning aid for the maxillary joint image**


Correct image acquisition is only possible with the chin support for maxillary joint images.

› Inserting the chin support for the maxillary joint image.


**Inserting the positioning aid for sinus images**

› Inserting the chin support for sinus images.

## Opening the head supports

- › Open the head supports by pressing the "Close/open head supports" button on the touch screen.
- › Open the cushion holder with cushions properly such that the patient can be positioned.

## 7.4 Positioning the patient

For the X-ray image, the patient is accurately aligned using the positioning beams.

Requirements:

- ✓ The patient has taken off jewelry and metal objects, e.g. earrings, hair slides, glasses, artificial dentures or orthodontic aids.
- ✓ The patient has put on a protective lead apron.
- ✓ The patient has been informed about the X-ray procedure.
- ✓ The patient has been informed that the unit may pass by close to his/her head (including through the field of vision). If patients feel uncomfortable with this, he or she can close his/her eyes while the image is being taken.
- ✓ The patient has been informed that he/she can press the emergency stop switch in the event of anxiety during image acquisition.
- ✓ The patient has been informed that he/she has to place his/her tongue against the roof of the mouth during the X-ray.
- ✓ The patient has been informed that he/she has to keep his/her eyes closed during positioning of the X-ray positioning beam.
- ✓ The patient has been told not to move while the X-ray is being taken until the unit is back in the starting position.

### CAUTION


**Danger of injury by moving rotating unit arm**

After the unit is switched on and the parameters are confirmed on the touch screen, the C-shaped angle connector piece is positioned. Persons can be injured during this.

- › Nobody may be present in the area of the rotating unit arm while the unit is being switched on.

- › Move the patient into an upright position at the unit.

It is also possible to position patients in a seated position (e.g. wheelchair users, tall patients).


- › Use the   buttons to set the height of the unit.

### CBCT image acquisition

The patient is positioned as follows depending on the indication:


- › The patient bites onto the bite block, with the upper and lower incisors resting in the grooves provided.

- › Use of the chin support for edentulous patients. Here, the patient places his/her chin on the chin support.

### Panoramic image acquisition


- › The patient bites onto the bite block, with the upper and lower incisors resting in the grooves provided.

- › Use the chin support for edentulous patients in the case of patients who do not have any teeth. Here, the patient places his/her chin on the chin support.


### Maxillary joint image

- › Position the patient with the upper lip against the chin support.

### Sinus image acquisition

- › Position the patient so that his/her bottom lip presses lightly against the chin support.

### Adjusting the position with the positioning beams

- i** The positioning beams use class 1 laser beams. Although these can dazzle the patient, they are safe and will not damage the eye.


- i** For CBCT images it is sufficient to perform positioning based on the mid-sagittal plane.

For all other image types the patient needs to be positioned more accurately with the aid of the following steps.

The alignment of the upper canine X-ray positioning beam is crucial for the image quality of panoramic images.

- › Check to make sure that the patient has closed his/her eyes.
- › If necessary, correct the height of the unit again.
- › Activate the positioning beams with the  button.

- › Check the X-ray positioning beam for the mid-sagittal plane and correct the position of the patient if necessary.

- › Align the head of the patient according to the Frankfurt horizontal plane with the aid of the X-ray positioning beam.

Exception: sinus image. Patient over-extends the cervical vertebral column to the rear by approx. 10° to 15°.




Fig. 3: Frankfurt plane: laser height at the lower edge of the eyes


- › Correct any inclination of the head via the height of the unit. If necessary, correct the X-ray positioning beam manually.

- › Direct the "upper canine" X-ray positioning beam as exactly as possible to the middle of the upper canine. Have the patient smile so the upper canine is visible.


- › If necessary, correct the X-ray positioning beam manually.

- › Once the patient has been correctly positioned using the positioning beams, deactivate the positioning beams using the  button.


**Adjusting the head supports**

- › Adjust the height of the head supports.

- › Carefully press the head supports by hand towards the head in order to check if they are in the right position. The device or the head supports are not damaged in the process. Ideally, the head supports should make contact slightly above the eye brows; correct the position as required.

- › Close the head supports with the  button. To do this, just press the button briefly – don't press and hold.

The head supports automatically touch against the head of the patient at a defined pressure.

**7.5 Start a test run**

The test run ensures that the unit can perform the image acquisition without any problems. This prevents unnecessary exposure of the patient to radiation.

-  No radiation is generated during the test run.

Requirements:

- ✓ The patient is correctly positioned in the unit using the positioning aids and the positioning beams.
- ✓ The required imaging program has been selected.
- › Touch **Test run** on the touch screen.
- › Tap **Run** and hold. While doing this, constantly monitor the movements of the unit. If the unit is obstructed in its movements, release the **Run** button. The unit stops immediately. Reposition the patient.
- › Press **Return run** to perform the return run.

## 7.6 Taking the X-ray image

### CAUTION

#### Injuries due to X-rays

X-rays can cause tissue damage.

- › Comply with the radiation protection regulations.
- › Maintain the minimum distance.


### CAUTION

#### Danger of excessively high radiation dose

- › Before an image acquisition is triggered, all data entered on the PC must be checked on the touch screen.

- › Check all parameters on the touch screen and change them if necessary. The changed parameters are immediately synchronized with the imaging software. The parameters can then no longer be changed in the imaging software.
- › Remind the patient to press his/her tongue against the roof of the mouth during image acquisition.

- › Press **Start** to confirm the parameters. The rotating unit is being positioned. The LED on the exposure switch and the status LED on the unit light up green. The touch screen displays that the unit is ready to take an image.


- › Trigger the image by pressing and holding the button on the exposure switch until the acoustic signal stops and the control lamp goes out. The scan times depend on the patient type, imaging program and image quality (see "12 Program parameters"). Image acquisition is started. During the acquisition, the LED on the exposure switch and on the unit illuminates yellow. An acoustic signal is issued.

If the button on the exposure switch is released before the control lamp goes out or the emergency stop switch is pressed (e.g. if there is a danger to the patient or to anyone else in the area), the ongoing image acquisition is aborted. The X-ray image will be unusable as a result and should be retaken as required. In this case, the operators must use their skills and training to decide on the risks of a repeated image acquisition.

In addition, an error messages appears on the touch screen.

While an X-ray is being taken, this is indicated on the touch screen by:


After the acquisition of maxillary joint images, it is necessary to confirm a message on the touch screen and trigger a second image acquisition. The images are then combined into a single image.

The LED on the unit lights up blue when the X-ray acquisition has been completed.

The rotating unit does not automatically move back to the starting position after the trigger button is released.

- › Click **OK** to confirm the message.

- › Release the head supports. The patient can leave the X-ray room.
- › Remove the hygienic protective cover.
- › Remove and disinfect the positioning aids.
- › The unit can be returned to its starting position by touching **Start position**. Otherwise, the C-shaped arm is positioned via the imaging software during adjustment of the parameters.

## 7.7 Emergency stop switch

The emergency stop switch stops the unit and switches it off. It can be used if the unit is taking an X-ray even though the exposure button is no longer being pressed, or if the patient is injured or

the unit is damaged. It can also be used to avert an unwanted collision.

The yellow stickers on the patient positioning system with the symbol  show the location of the emergency stop switch.


Fig. 4: Emergency stop switch sticker on the operator side



Fig. 5: Emergency stop switch sticker on the patient side


Misuse of the emergency stop switch can lead to data loss.

- › Press the emergency stop switch.

Device is switched off.

The unit automatically restarts.

### Releasing the emergency stop switch

#### CAUTION

#### Danger of injury by moving rotating unit arm

After the unit is switched on and the parameters are confirmed on the touch screen, the C-shaped angle connector piece is positioned. Persons can be injured during this.

- › Nobody may be present in the area of the rotating unit arm while the unit is being switched on.

- › Twist the emergency stop switch to release it.

## 8 Cleaning and disinfection

### NOTICE

**The use of unsuitable agents and methods can damage the unit and accessories**

- › Only use the disinfection and cleaning agents specified or approved by Air Techniques and the EPA.
- › Comply with the operating instructions of the disinfectants and cleaning agents.

Wear hand protection.

Prior to working on the unit or in case of danger, disconnect it from the mains.

The following cleansers and disinfection agents have been tested for the compatibility of materials:

- Birex wipes
- Discide Ultra Towelettes
- Opti Cide 3 surface wipes
- Optim 33TB wipes
- Maxiwipe germicidal cloth
- Cavi Wipes
- Cavi Wipes 1
- PDI Sani-Cloth
- Clorox
- Monarch

### 8.1 Surface of the unit

### NOTICE

**Damage to the touch screen caused by cleaning it with disinfectant**

- › Only clean the touch screen with a soft cloth and a commercially available cleaning agent.

The surface of the unit must be cleaned and disinfected if it is contaminated or soiled.

### NOTICE

**Liquid can cause damage to the unit**

- › Do not spray the unit with cleaning agents or disinfectants.
- › Make sure that liquid penetrates into the unit.


- › Remove any soiling with a soft, damp, lint-free cloth.
- › Disinfect the surfaces with a disinfection wipe. Alternatively, use a quick-acting surface disinfectant on a soft, lint-free cloth. Comply with the operating instructions of the disinfectant.

### 8.2 Positioning aids


The positioning aids must be cleaned and disinfected if they are contaminated or soiled.

#### Head supports with cushion

- › Pull the head supports off the device.
- › Remove the cushions from the head supports.

- › Remove the cushion holder.

- › Remove any soiling with a soft, damp, lint-free cloth.

- › Disinfect the surfaces using a disinfection wipe. Alternatively, use a quick-acting surface disinfectant on a soft, lint-free cloth. Comply with the operating instructions of the disinfectant.
- › Reprocessing the cushions (see "9 Reprocessing").

## 9 Reprocessing

The following accessories need to be reprocessed:

- Bite block:
  - Manual cleaning
  - Automatic cleaning and disinfection
  - Steam sterilization
- Adapter bite block, chin holder for mandibular joint image, chin holder for edentulous jaws, and chin holder for sinus image
  - Manual cleaning
  - Manual disinfection
  - Automatic cleaning and disinfection
- Cushion for head supports Plus
  - Manual cleaning
  - Manual disinfection
  - Automatic cleaning and disinfection

In order to prevent damage to the accessories, only the methods described above must be used.

### 9.1 Risk analysis and classification

A risk analysis and classification of medical devices that are common in dentistry must be performed before they are reprocessed by the operator. Comply with the country-specific guidelines, standards and regulations.

Accessories of the medical device are also subject to reprocessing.

#### Recommended classification according to EN ISO 17664, Appendix C

Recommended classification given proper use of the bite block:

##### **semi-critical**

Recommended classification given proper use for the adapter for bite block, chin holder for mandibular joint image, chin holder for edentulous jaws, and chin holder for sinus image and cushions for head supports Plus:

##### **non-critical**

The operator is responsible for correct classification of the medical products, defining the reprocessing steps and performing the reprocessing.

### 9.2 Reprocessing procedure in accordance with EN ISO 17664

Perform the reprocessing procedure after each patient treatment and according to the reprocessing procedure established by EN ISO 17664.

#### **Important information!**

The reprocessing instructions in accordance with EN ISO 17664 have been independently tested by the manufacturer for the preparation of the device and its components for their reuse.

The person conducting the reprocessing is responsible for ensuring that the reprocessing is performed using equipment, materials and personnel that attains the desired results. This requires validation and routine monitoring of the reprocessing process. Any deviation from the instructions described herein by the staff reprocessing the equipment could lead to lower effectiveness and possible negative consequences: these lie solely with the staff responsible.

Frequent reprocessing has little effect on the components of the device. The end of the product life cycle is mainly influenced by the amount of wear and tear or damage resulting from its use.

The use of soiled, contaminated and damaged components is at the sole responsibility of the person performing the reprocessing and the operator.

The reprocessing procedure was validated as follows:

- Pre-cleaning:
  - Monarch surface disinfection wipes (Air Techniques)
  - Cleaning brush
- Manual cleaning:
  - Monarch enzymatic cleaner 2.0 % (v/v) (Air Techniques)
- Manual disinfection:
  - Monarch surface disinfection wipes (Air Techniques)
- Automated cleaning and disinfection was performed in accordance with EN ISO 15883 with tested efficacy:
  - Washer-disinfector PG 8536 CD (Miele, Gütersloh, Germany)
  - Cleaning agent: Neodisher MediClean Forte
  - Programs: *Cleaning without neutralization* and *THERMAL DES*
- Steam sterilization:
  - Steam sterilizer Systec DX-45 (Systec GmbH, Linden, Germany)

### 9.3 General information

- › Comply with all national directives, standards and specifications for the cleaning, disinfection and sterilization of medical devices as well as the specific specifications for dental practices and clinics.
- › Comply with the specifications (in "9.5 Manual cleaning, intermediate rinsing, disinfection, drying" and "9.6 Automatic cleaning, intermediate rinsing, disinfection, final rinse, drying") when you select the cleaning agents and disinfectants to be used.
- › Comply with the concentrations, temperatures, residence times and post-rinsing specifications issued by the manufacturer of the cleaning agent and disinfectant.
- › Only use cleaning agents that are non-fixing and aldehyde-free and display material compatibility with the product.
- › Only use disinfectants that are aldehyde-free and display material compatibility with the product.
- › Do not use any rinse aid (danger of toxic residue on the components).
- › Only use freshly-produced solutions.

- › Only use distilled or de-ionized water with a low bacterial count (at least drinking water quality) that is free from facultatively pathogenic micro-organisms (e.g. Legionella bacteria).
- › Use clean, dry, oil- and particle-free compressed air.
- › Do not exceed temperatures of 281 °F.
- › Subject all devices used (e.g. ultrasonic bath, cleaning and disinfection device (washer-disinfector), sealing device, steam sterilizer) to regular maintenance and inspections.

### 9.4 Preparation at the operating location

Wear hand protection.

Wear eye protection.

Use a mask.

Wear protective clothing.

#### WARNING

#### Risk of infection from contaminated products

Risk of cross contamination

- › Reprocess the product correctly and promptly before its first use and after every subsequent use.
- › Protect the unit from contamination when transporting it from the treatment chair to the reprocessing location.
- › Brush off all surfaces below the water surface with a soft hygienic brush until they are clean to the eye.
- › Wipe off all surfaces for at least one minute with a disinfection wipe.

### 9.5 Manual cleaning, intermediate rinsing, disinfection, drying

#### Cleaning

- › Place individual parts in a cleaning agent bath making sure that all parts are covered.

- › Comply with the exposure times of the cleaning agent.

### Intermediate rinsing

After the action time prescribed by the manufacturer has expired:

- › Rinse all components with water for at least 1 minute (temperature < 95°F).

### Disinfecting

- › Wipe off all surfaces for at least one minute with a disinfection wipe.

### Drying

- › If necessary, re-dry at a clean location using a hygienic, lint-free cloth.
- › Blow dry the components with compressed air in a clean location.

## 9.6 Automatic cleaning, intermediate rinsing, disinfection, final rinse, drying

### Selection of the washer-disinfector

Automatic cleaning and disinfection requires a washer-disinfector with the following properties and validated processes:

- Corresponds to and tested in accordance with EN ISO 15883
- Certified program for thermal disinfection ( $A_0$  value  $\geq 3000$  or at least 5 minutes at 382°F)

Program is suitable for the components and includes sufficient rinsing cycles.

For more information: "9.3 General information".

### Selection of the cleaning agents for automatic cleaning

The following properties are required:

- Material compatibility with the product
- Compliance with the washer-disinfector manufacturer's specifications

For further information, see "9.3 General information".

### Cleaning and disinfection

- › Place all components in the washer-disinfector (follow the manufacturer's instructions).
- › Make sure there are no hidden areas that are missed by the rinsing process.
- › Secure the components with a suitable fixture of the washer-disinfector.

## 9.7 Check for function

- › After the end of the cleaning and disinfection cycle, check the components for any residual soiling and residual moisture. If necessary, repeat the cycle.
- › If necessary, replace any damaged parts.
- › The components should be packaged as soon as possible after drying and checking.

## 9.8 Steam sterilization

### Packing

For packing of the components, only use transparent paper-film sterilization packaging that is approved for use in steam sterilization according to the manufacturer's instructions. This includes:

- Temperature resistance up to 281°F
- Standards DIN EN ISO 11607-1/2
- The applicable sections of the standard series DIN EN 868

The sterilization packaging must be sufficiently large. Once it is loaded, the sterilization packaging may not be under any strain.

### Steam sterilization

#### WARNING

#### Incorrect sterilization reduces effectiveness and can damage the product

- › Only steam sterilization is permitted.
- › Comply with the specified process parameters.
- › Comply with the manufacturer's instructions regarding use of the steam sterilizer.
- › Do not use any other methods.

### Steam sterilizer requirements:

- Complies with EN 13060 or EN 285 and/or ANSI AAMI ST79
- Suitable programs for the products listed (e. g. with hollow bodies: fractionated vacuum procedure including three vacuum steps)
- Sufficient drying of the product
- Validated processes in accordance with ISO 17665 (valid IQ/OQ and product-specific performance appraisal (PQ))

Perform the following steps:

- › Sterilize the parts for sterilization (at least 20 minutes at 250°F, at least 4 minutes at 270°F or at least 5 minutes at 274°F).

 Do not exceed 281 °F.

### Marking

- › Mark the packaged, treated medical device appropriately such as to ensure its safe application.

## 9.9 Issue clearance for the parts for sterilization

The reprocessing of the medical device ends with the documented clearance for storage and renewed use.

- › Document the release of the medical device after reprocessing.

## 9.10 Storing parts for sterilization

- › Comply with the stated storage conditions:
  - Store the parts protected against contamination
  - Dust-protected, e.g. in a locked cabinet
  - Protected against moisture
  - Protected against excessive temperature fluctuations
  - Protected against damage

The integrity of the packaging of a sterile medical device be lost as a result of a particular incident and the passage of time.

Potential external contamination of the sterile barrier system should be taken into account in terms of aseptic preparation when establishing the storage conditions.

# 10 Maintenance

## 10.1 Recommended maintenance schedule

Only qualified personnel or personnel trained by Air Techniques is allowed to service the unit.

**DANGER**  
**Risk of electric shock due to live parts**

- › Install an all-pole disconnect switch with a contact opening width of at least 3 mm in the electrical connection to the mains power supply. It must be possible to secure the disconnect switch so that it cannot be inadvertently switched back on again.

**WARNING**  
**Risk of infection from contaminated products**

Risk of cross contamination

- › Reprocess the product correctly and promptly before its first use and after every subsequent use.

**NOTICE**  
**Damage to the X-ray tube due to overheating**

- › Observe the cooling curves of the X-ray tube when working with the service tool.

Inspection interval	Inspection work
Every 3 years	<ul style="list-style-type: none"> <li>› Functional test of the display. Are all symbols displayed?</li> <li>› Functional test of the exposure button.</li> <li>› Do the various status LEDs light up?</li> <li>› Check that the mechanism of the head supports works correctly. Are the head supports easy to detach and put back on?</li> <li>› Functional test of the EMERGENCY OFF button. Is the EMERGENCY OFF button easy to operate mechanically, and does it light up when pressed?</li> <li>› Light barrier test for all light barriers installed in the unit.</li> <li>› Visually check the positioning beams. Check the operation of the levers for adjusting the beam localizers.</li> <li>› Check the X-ray images for artifacts. If necessary, adjust the collimator and/or calibrate the sensor.</li> <li>› Check the firmware and software versions.</li> <li>› Perform a comparative dose measurement based on the requirements from the acceptance test.</li> <li>› Recurring tests and tests after repairs to medical electrical equipment – IEC 62353 (VDE 0751-1).</li> </ul>

Maintenance interval	Maintenance work
Every 3 years	<ul style="list-style-type: none"><li data-bbox="297 169 972 193">› Visually and acoustically check the linear movement on the rotating unit.</li><li data-bbox="297 204 1021 252">› Check the operation of the lift motor. Does the unit lift and lower with minimal noise?</li></ul>

# ? Troubleshooting

**CAUTION**

Any oil leaking from the X-ray tube in the event of a fault is harmful.

- › Wipe up any oil immediately.
- › Do not swallow the oil.
- › Stop using the unit and inform a service technician.

## 11 Tips for operators and service technicians

Any repairs above and beyond routine maintenance may only be done by suitably qualified personnel or by one of our service technicians.

Error	Possible cause	Remedy
<b>Unit does not switch on</b>	EMERGENCY STOP SWITCH accidentally activated	› Releasing the EMERGENCY STOP SWITCH.
	No mains voltage	› Check the mains cable and electrical connection; replace if necessary. › Contact technician.
		› Check the mains fuse in the building.
	On/off switch is defective	› Contact technician.
<b>Unit not responding</b>	The unit has not yet completed the startup procedure	› After switching on, wait until the booting process has finished.
	Cables not correctly connected	› Check the cable connections.
	Plug-in contacts of the fiber optic cable are contaminated	› Clean plug-in contacts and sockets.
	Driver for PCI Express frame grabber card is not installed or incorrectly installed	› Install driver or complete VistaVox Plugin once again.
	COM port incorrectly configured	› Check the COM port in the service tool.

Error	Possible cause	Remedy
<b>Error messages during start of the acquisition of an X-ray image or during shut down of the PC</b>	Energy-saving options incorrectly configured	› Deactivate the energy-saving options in Windows and the BIOS completely.
	Supply voltage for graphic card inadequate or incorrectly connected	› Check the plug connections. › Compare requirements of the graphic card with the power supply of the PC, use larger power supply if necessary.
	PC and/or graphic card fail to comply with the specified system requirements	› Set up the system in accordance with the system requirements.
	User account control (UAC) has not been correctly configured	› Make the user account control settings according to the information in the installation instructions.
	USB dongle not detected	› Check whether the USB dongle (included in scope of delivery) is plugged into the reconstruction computer, or check that the USB dongle is correctly plugged in.
	Virus scanner prevents the acquisition of an X-ray image	› Add the installation paths of the imaging software as exceptions in the virus scanner.
	Firmware of the device does not match the software version	› Check software versions and update if necessary.
	The device calibration has not been imported or incompletely imported	› Carry out initial commissioning using the service tool again/initially.
Door contact is not closed	› Check door contact and plug-in connections of the door contact, close door properly.	

## 12 Program parameters

The digital extraoral dental X-ray system meets the requirements set out in standard IEC 60601-2-63. The dosage information complies with the requirements of the standard and is stated in mGy.

The accuracy of the DAP/dose values is  $\pm 50\%$ .

### 12.1 CBCT program parameters

CBCT image acquisition, normal image volume, 16.4 s								
	4.0 mA		6.3 mA		8.0 mA		10.0 mA	
	mGy	mGycm <sup>2</sup>	mGy	mGycm <sup>2</sup>	mGy	mGycm <sup>2</sup>	mGy	mGycm <sup>2</sup>
75 kV	4.00	215.92	6.13	331.05	7.79	420.70	9.74	526.41
79 kV	4.50	242.90	6.89	372.41	8.76	473.26	10.96	592.17
90 kV	6.00	324.11	9.20	496.93	11.69	631.50	14.63	790.18
94 kV	6.55	353.92	10.04	542.63	12.77	689.58	15.97	862.85

CBCT image acquisition, normal image volume 5x5, 11 s										
	4.0 mA		6.3 mA		8.0 mA		10.0 mA		11.0 mA	
	mGy	mGyc m <sup>2</sup>	mGy	mGyc m <sup>2</sup>	mGy	mGyc m <sup>2</sup>	mGy	mGyc m <sup>2</sup>	mGy	mGyc m <sup>2</sup>
79 kV	3.02	118.85	4.63	182.22	5.88	231.56	7.36	289.75	8.08	318.22
94 kV	4.40	173.17	6.74	265.51	8.57	337.41	10.72	422.19	11.78	463.68
98 kV	4.77	187.76	7.31	287.87	9.29	365.83	11.63	457.75	12.77	502.73

### 12.2 Panoramic program parameters

Pano image acquisition, normal jaw arch, normal patient, quality HQ, 13.5 s						
	4.0 mA		6.3 mA		10.0 mA	
	mGy	mGycm <sup>2</sup>	mGy	mGycm <sup>2</sup>	mGy	mGycm <sup>2</sup>
60 kV	3.85	24.28	6.06	38.15	9.57	60.31
67 kV	4.74	29.84	7.43	46.83	11.76	74.09
70 kV	5.12	32.24	8.03	50.59	12.70	80.03
74 kV	5.66	35.66	8.88	55.95	14.05	88.52
80 kV	6.47	40.79	10.16	64.00	16.07	101.25

	12.5 mA		14.0 mA	
	mGy	mGycm <sup>2</sup>	mGy	mGycm <sup>2</sup>
60 kV	11.77	74.18	13.21	83.23
67 kV	14.46	91.07	16.22	102.21
70 kV	15.61	98.37	17.52	110.40

	12.5 mA		14.0 mA	
	mGy	mGycm <sup>2</sup>	mGy	mGycm <sup>2</sup>
<b>74 kV</b>	17.27	108.81	19.38	122.11
<b>80 kV</b>	19.76	124.46	22.17	139.68

## 13 Information on scattered radiation

### 13.1 CBCT scattered radiation

Test equipment: Dosimeter Radcal 9015

Test conditions	
Program parameters	CBCT
Image acquisition volume	Normal
Voltage	99 kVp
Current	14 mA

R °	1 m mR/h	1.5 m mR/h	2 m mR/h
0	588.2	135.3	87.1
45	549.3	249.4	106.8
90	472.6	307.3	78.4
135	458.8	287.6	89.3
180	12.9	4.6	1.3
225	410.5	288.7	98.2
270	663.2	301.4	112.4
315	429.7	194.2	92.3

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
For over 50 years, Air Techniques has been a leading innovator and manufacturer of dental products. Our priority is ensuring complete satisfaction by manufacturing reliable products and providing excellent customer and technical support. Whether the need is digital imaging, utility room equipment or merchandise, Air Techniques can provide the solution via our network of authorized professional dealers.

Air Techniques' family of quality products for the dental professional include:

- **Digital Imaging**
  - Digital Radiography
  - Intraoral Cameras
  - Caries Detection Aid
  - X-ray Systems
  - Film Processors
- **Utility Room**
  - Dry Vacuums
  - Wet Vacuums
  - Air Compressors
  - Amalgam Separator
  - Utility Accessories
  - Utility Packages
- **Merchandise**
  - Imaging Accessories
  - Chemistry
  - Processor Cleaners
  - Surface Disinfectant
  - Instrument Cleaner
  - Hand Sanitizer + Hand Lotion
  - Evacuation System Cleaner
  - Water Line Cleaner

**Manufactured for / Distributed by:**

Air Techniques, Inc. | 1295 Walt Whitman Road | Melville, New York 11747-3062, USA  
Phone: 800-247-8324 | Fax: 888-247-8481

DÜRR DENTAL SE | Höpfigheimer Str. 17 | 74321 Bietigheim-Bissingen | Germany

[www.airtechniques.com](http://www.airtechniques.com)

**AIR**  
**TECHNIQUES** equipped for life®