


CLINICAL AIR COMPRESSOR MAINTENANCE KIT

Please contact your preferred authorized dealer for purchase (PN E5256)



Note: - - - - -
Indicates part is located within device

(Two additional Relief Valve Mufflers and Muffler Intakes are located on the opposite side)







Replace every year **Replace every 2 years**

1	2	3	4	5	6
E5253 Fan Filter Qty: 2/year	E5131 Bacteria Filter Qty: 1	E5319 Relief Valve Muffler Qty: 4	C1162 Muffler Intake Qty: 4	E5207 Battery Qty: 1	E5254 Blue Relay Qty: 5
					



Click Here or scan QR code to access videos detailing parts replacements

REPLACEMENT PART TRACKER FOR CLINICAL COMPRESSOR SYSTEMS

Air Compressor System Serial Number: _____						
REPLACE EVERY YEAR						
Image	Item	PN	Date	Date	Date	Date
	1	E5253				
REPLACE EVERY 2 YEARS						
	2	E5131				
	3	E5319				
	4	C1162				
	5	E5207				
	6	E5254				
Signature: _____						

Clinical Installations from Leading Dental Facilities around the World



We are here to support the seamless integration of your entire project from the planning phase through after sales training.

While your preferred authorized dealer will provide the service you require, please don't hesitate to contact us directly for additional customer service and technical support.

phone: 800.247.8324 opt 2 for tech support.

email: techservice@airtechniques.com

www.airtechniques.com

Corporate Headquarters

1295 Walt Whitman Road | Melville, NY 11747-3062 | Phone 800.247.8324 | Fax 888.247.8481



CORPORATE HEADQUARTERS
1295 Walt Whitman Road
Melville, New York 11747- 3062
800-247-8324
www.airtechniques.com

Air Techniques, Inc.
© Copyright 2024
PN E5256-LIT Rev. B

