


CLINICAL VACUUM MAINTENANCE KIT

Please contact your preferred authorized dealer for purchase (PN E5257)



Note: ----- indicates part is located within device



Replace when required*

Replace every year or when full†

Replace every year

C

57630
CleanStream
Qty: 1



1

E5148
Cassette
Qty: 2/year



2

E5253
Fan Filter
Qty: 2/year




3

9000320002
Gasket
Qty: 2/year



4

E5133
Non-Return Flap CST
Qty: 1/year



Notes: (* & †)-see the following page

Replace every 2 years

5

E5142
Tank Float Switch
Qty: 3



6

E5286
Cleaner Float Switch
Qty: 1



7

E5207
Battery
Qty: 1




8

9000341004
Drain Pump Impeller
Qty: 1



9

54139
Pump Check Valve
Qty: 6




10

E5136
Coarse Filter CST
Qty: 1



11

E5255
Orange Relay
Qty: 1



12

E5254
Blue Relay
Qty: 1

















13

E5131
Bacteria Filter
Qty: 2




Click Here or scan QR code to access videos detailing parts replacements

REPLACEMENT PART TRACKER FOR CLINICAL VACUUM SYSTEMS

Vacuum System Serial Number:										
Image	Item	PN	Date	Date	Date	Date				
REPLACE WHEN REQUIRED										
	C	57630*								
REPLACE EVERY YEAR OR WHEN FULL										
	1	E5148†								
REPLACE EVERY YEAR										
	2	E5253								
	3	9000320002								
	4	E5133								
REPLACE EVERY 2 YEARS										
	5	E5142	X	X	X	X				
	6	E5286								
	7	E5207								
	8	9000341004								
	9	54139								
	10	E5136								
	11	E5255								
	12	E5254								
	13	E5131								
Signature:										

Notes:

†Amalgam Cassette (PN E5148): To stay in compliance with EPA regulations, go to (<https://www.airtechniques.com/login/>) to register your amalgam separator, generate the shipping label and retrieve your recycling certificate. (For USA and Canada only).

*** Item purchased separately:**

 CleanStream (PN 57630)-Monarch CleanStream: The only cleaner we recommend to ensure the highest performance and longevity of your equipment.

Clinical Installations from Leading Dental Facilities around the World



We are here to support the seamless integration of your entire project from the planning phase through after sales training.

While your preferred authorized dealer will provide the service you require, please don't hesitate to contact us directly for additional customer service and technical support.

phone: 800.247.8324 opt 2 for tech support.
email: techservice@airtechniques.com

www.airtechniques.com

Corporate Headquarters

1295 Walt Whitman Road | Melville, NY 11747-3062 | Phone 800.247.8324 | Fax 888.247.8481



CORPORATE HEADQUARTERS
 1295 Walt Whitman Road
 Melville, New York 11747- 3062
 800-247-8324
www.airtechniques.com

Air Techniques, Inc.
 © Copyright 2024
 PN E5257-LIT Rev. B

